


# Para-dox<sup>®</sup>

Photoredox Catalysis &  
Parallel Synthesis Reaction Blocks



## Parallel Synthesis

## Photoredox Catalysis



## Electrochemistry

Winter 2023-24



analytical

# Table of Contents

## Parallel Synthesis

Aluminum Reaction Blocks.....2

Standard (1mL) 24-Well: 2 • Standard (1mL) 96-Well: 3  
HPLC (2mL) 48-Well: 4 • 1 DRAM (4mL) 24-Well: 5 • 2 DRAM (8mL), 24-Well: 5



2

## NEW! Gen II Reactors



## Photoredox Catalysis

NEW! Gen II Aluminum Reaction Blocks.....6

Aluminum Reaction Blocks.....8

Micro (50µL) 24-Well: 8 • Micro (50µL) 96-Well: 9 • Standard (1mL) 24-Well: 10 • Standard (1mL) 96-Well: 11  
HPLC (2mL) 48-Well: 12 • 1 DRAM (4mL) 24-Well: 13 • 2 DRAM (8mL), 24-Well: 13



Temperature Controlled Reactors (TCR).....16

Lumidox Gen II LED arrays for Photoredox.....18

Controller: 19 • 96-position Arrays: 20 • Discovery: 23 • Other Arrays: 24

Thermal Transfer Decks (TTD).....25

Cell Culture / PCR Adapters .....26

Flow Reactor for Photoredox.....28

LumLamp for Photoredox.....29



6

## Electrochemistry

HTe-Chem .....30

Flow Electrolysis.....32



## Accessories

EquaVAP® Evaporators: 34 • Collection Plates and Cap Mats: 36 • Thermal Adapter Plates: 37  
Filter Plates and Vacuum Manifold: 38 • Vial Trays and Loaders: 39 • Powder transfer Plates: 40



30

34



Photoredox Catalysis &  
Parallel Synthesis Reaction Blocks

**Para-dox® Aluminum Reaction Blocks**

A successful approach to high-throughput reaction screening requires the best tools available. Our comprehensive line of Photoredox Catalysis and Parallel Synthesis reaction blocks are exactly what you need. Designed specifically for high-throughput screening (HTS) applications, the SBS (SLAS/ANSI) format allows for use in any industry standard automation application.

- Useful for generating compound arrays (medicinal chemistry)
- Useful for conducting screening reaction conditions
- Can be used on tumble stirrers / hot plate stirrers / robotic platforms
- Can easily be used in a glovebox


- Temperature range -78°C – 150°C
- Validated to have less than 5% solvent loss with prolonged heating above boiling point
- Silicone Rubber Mats provide compression sealing
- Teflon® PFA Films keep the glass reaction vials from sticking to the silicone rubber mats during heating

*Tech Note: Parallel Synthesis Reaction Blocks are ideal for the optimization of chemical processes, lead generation optimization, and screening for optimal reaction conditions.*

**Parallel Synthesis Reaction Blocks  
for Optimization of Chemical Processes**




24249



24240

**Standard 24-Position Parallel Synthesis Reaction Block  
for 8x30mm, 1mL Glass Inserts**



| Cat. No. | Description   | Qty  |
|----------|---|------|
| 24240    | Parallel Synthesis 24-Well Block Assembly. Includes Vial Rack, Top Cover, Bottom Cover, PFA Films, Rubber Mats, Screws and Rack Adapter | Each |
| 24249    | Parallel Synthesis/Optimization 24-Well Block Assembly, <i>no Rack Adapter</i>  | Each |
| 24245    | Rack Adapter  | Each |
| 24247    | Replacement Films for 24 Position, Optimization Block   | 25   |
| 24248    | Replacement Mats for 24 Position, Optimization Block  | 25   |
| 3/8\"/>  |   |      |

**Accessories for  
24-Well Parallel Synthesis Blocks**

| Cat. No.   | Description  | Qty  |
|------------|--|------|
| 84001-Case | 1mL Clear Glass Shell Vials, 8 x 30mm  | 1000 |
| 13258      | Stainless Steel Cylinder Stir Bars, 1.98mm x 4.80mm (for 8x30mm Vials)   | 1000 |
| SD1000     | Milwaukee M4 1/4in. Hex Screwdriver Kit. Includes 2 Batteries and Charger (Note: batteries not available for international shipping) | Each |



SD1000



96960

| Cat. No. | Description  | Qty  |
|----------|--|------|
| 96960    | Parallel Synthesis 96-Well Block Assembly Includes Vial Rack, Top and Bottom Covers, PFA Films, Rubber Mats and Screws | Each |
| 96967    | Replacement Films for 96-Well Optimization Block   | 25   |
| 96965    | Replacement Mats for 96-Well Optimization Block  | 25   |
| 1/2\"/>  |  |      |

**Recommended for top performance!  
Assembled Vials in Stackable Trays**

- For easy loading - *Saves Time!*
- Our thorough QC process ensures *Less Evaporation* compared to loose vials



884001

(with 8x30mm Shell Vials)

| Cat. No. | Description   | Qty  |
|----------|---|------|
| 884001   | Assembled Stackable Tray Loaded with 8x30mm Shell Vials (84001-CASE), Includes Tray and Vials | Each |




96242

(with 5x31mm Flat-Bottom Vials. Must use with Aluminum Spacers, part # 96969)

| Cat. No. | Description   | Qty  |
|----------|---|------|
| 96242    | Assembled Stackable Tray Loaded with 5x31mm Flat-Bottom Vials (20303-CASE), Includes Tray and Vials | Each |

**Standard 96-Position Parallel Synthesis  
Reaction Block  
for 8x30mm, 1mL Glass Inserts**


Replace PFA Films with every use for best performance



1\"/>

1/2\"/>

Both the 24-well and 96-well blocks will hold 750µL, 8x30mm flat-bottom glass inserts (part no. 84001). They will also hold 250µL, 5x31mm inserts (part no. 20303) when used with anodized aluminum spacers (part no. 96969).




20303

84001

96969

**Accessories**

| Cat. No. | Description  | Qty  |
|----------|--|------|
| 884001   | Assembled Tray with 8 x 30mm Shell Vials   | Each |
| 96969    | Anodized Aluminum Spacers for 5 x 31mm Inserts   | 96   |
| 13258    | SS Cylinder Stir Bars, 1.98mm x 4.80mm (for 8x30mm Vials)  | 1000 |
| SD1000   | Milwaukee M4 1/4in. Hex Screwdriver Kit, includes 2 Batteries and Charger (Note: batteries not available for international shipping) | Each |




**Top Cover with Wider Holes  
for 96-Well Reaction Blocks**

- Enlarged (4.3mm) Holes
- 316 Stainless Steel
- Made for Automated Sampling Instruments

| Cat. No. | Description   | Qty  |
|----------|---|------|
| 96985    | Reaction Block Cover with Enlarged Holes, 6mm Thickness | Each |
| 96986    | Reaction Block Cover with Enlarged Holes, 3mm Thickness | Each |


**HPLC 48-Position Parallel Synthesis Reaction Block**  
for 12x32mm, 2mL Glass Vials



Use with 12 x 32mm vials, sold separately.

| Cat. No.    | Description  | Qty  |
|-------------|--|------|
| 48012       | 48-Well Aluminum Reaction Block for 12mm (2mL) OD Vials. Includes: Base Plate, Cover, Screws, PFA Film and 2 Rubber Mats | Each |
| 48482       | Rubber Mats for 48 Well (12mm OD) Reaction Blocks  | 25   |
| 48483       | PFA films for 48 Well (12mm OD) Reaction Blocks  | 25   |
| VSCREW48    | 1 1/4" Screws for 48-Well Aluminum Reaction Plate  | 100  |
| 11211-Case  | 11mm, 12 x 32, 2mL Wide Mouth Glass Crimp/Snap Vials   | 1000 |
| 11221-Case  | 11mm, 12 x 32, 2mL Wide Mouth Glass Crimp Vials  | 1000 |
| 11221M-Case | 11mm, 12 x 32, 2mL Wide Mouth Glass Crimp Vials with Marking Spot  | 1000 |

**1 Dram, 24-Position Parallel Synthesis Reaction Block with 20mm Well Spacing**  
for 14x45mm, 4mL Glass Vials - fits on Mettlers QX96




Use with 15 x 45mm vials, sold separately. See part number 31531-Case.

| Cat. No. | Description   | Qty  |
|----------|---|------|
| 24015    | 24-Well Aluminum Reaction Block for 15mm OD (1 Dram) Vials, Includes: Base Plate, Cover, Screws, PFA Film and 2 Rubber Mats | Each |
| 24282    | High-Temp Silicone/Rubber Cover for 24-Well Reaction Blocks for Vials   | 25   |
| 24283    | PFA Film (0.005" Thick) for 24-Well Reaction Blocks for Vials   | 25   |
| VSCREW24 | 1 3/4" Screws for 24-Well Aluminum Reaction Plate   | 100  |

| Cat. No.   | Description   | Qty  |
|------------|---|------|
| 31531-Case | Advantage™ 13mm, 15 x 45, 4mL Clear Glass Screw Vials, Case | 1000 |
| 31554-Case | PTFE Lined Cap  | 1000 |




**SS Cylinder Stir Bars for 2mL HPLC Vials**

- Parylene Encapsulated
- For Photoredox and Parallel Synthesis Reaction Block Systems

| Cat. No. | Description   | Qty  |
|----------|---|------|
| 13258    | Stainless Steel Cylinder Stir Bars, 1.98mm x 4.80mm, for 2mL HPLC Vials | 1000 |

**Magnetic Stir Bars**

| Cat. No. | Description                                | Qty |
|----------|--|-----|
| 50225    | Magnetic Stir Bars, 5mm x 2mm, PTFE Coated | 25  |
| 502100   | Magnetic Stir Bars, 5mm x 2mm, PTFE Coated | 100 |



**48-Position Numbered Crimp Vial Assembly for Reaction Blocks**

| Cat. No. | Description  | Qty  |
|----------|--|------|
| 48221M   | 48-Well Numbered Crimp Vial Assembly. Includes: Vials (11221M) Marked 1-48; Stir Bar (13258) in Each Vial; 48-Well Cap Mat (99948) | Each |



**Cap Mat for HPLC (2mL) Vials in 48-Position Block**

| Cat. No. | Description  | Qty |
|----------|--|-----|
| 99948    | Clear Silicone/PTFE Cap Mat for HPLC Vials in 48-Position Reaction Block | 5   |


**Parallel Synthesis vs. Photoredox Catalysis Reaction Blocks**

Reactors for both applications are nearly identical. The difference lies in the base (vial rack) of the block. Reactors for parallel synthesis have wells with closed bottoms, whereas photoredox blocks have wells with open holes in the bottom allowing for light transmission via a Lumidox® II Array, as well as to read 2D barcodes.

Parallel Synthesis Block




Photoredox Catalysis Block



**NEW!** Para-dox® Gen II Reaction Blocks for Parallel Synthesis are now available. See website.

**2 Dram, 24-Position Parallel Synthesis Aluminum Reaction Block**  
for 17x60mm, 8mL Glass Vials



Use with 17 x 60mm vials, sold separately. See part number 31760-Case.

| Cat. No. | Description   | Qty  |
|----------|---|------|
| 24017    | 24-Well Aluminum Reaction Block for 17mm OD (2 Dram) Vials, Includes: Base Plate, Cover, Screws, PFA Film and 2 Rubber Mats | Each |
| 24282    | High-Temp Silicone/Rubber Cover for 24-Well Reaction Blocks for Vials   | 25   |
| 24283    | PFA Film (0.005" Thick) for 24-Well Reaction Blocks for Vials   | 25   |
| VSCREW24 | 1 3/4" Screws for 24-Well Aluminum Reaction Plate   | 100  |

| Cat. No.   | Description  | Qty  |
|------------|--|------|
| 31760-Case | Advantage™ 17 x 60mm, 8mL Clear Glass Screw Vials, Case                  | 1000 |
| 31542-CASE | 15mm Solid Black Polypropylene Screw Caps with PTFE/F217 Liners          | 1000 |
| 31543-CASE | Black Open Top PP Cap with PTFE/Silicone (0.065") Septa, 15-425mm Thread | 1000 |

**NEW! Gen II Reaction Blocks**

Para-dox SLAS footprint aluminum reaction blocks are now available in a 4 or 5-bolt configuration. Some benefits of these new blocks include:

- Less bolts = less time to assemble
- Less light bleed (48 and 96 position blocks)
- Less Maintenance
- Vials fit tighter (less shifting)
- Port for thermocouple, sensors, etc.
- Corner holes on both the lid and the bottom for accessory attachment
- Longer lasting screws (larger diameter, larger threads)
- Larger holes in lid to accommodate a variety of needle sizes
- Only one 1/8" rubber mat compared to two 1/16" mats on previous versions (makes for easier assembly)



**Photoredox Catalysis & Parallel Synthesis Reaction Blocks**



101960

**Standard 96-Position Photoredox Reaction Block, Gen II for 1mL, 8x30mm Glass Inserts**

| Catalog No. | Description   | Qty  |
|-------------|---|------|
| 101960      | 4-Bolt, 96-Well Aluminum Reaction Block for 8x30mm (1mL) Glass Inserts. Includes: Base Plate, Cover, Screws, PFA Film and Rubber Mats | Each |
| 101968      | Top PFA Films for Gen II 96-Well (101960) and 48-Well (101480) Reactors   | 25   |
| 101967      | Top Rubber Mats for Gen II 96-Well (101960) and 48-Well (101480) Reactors   | 25   |
| 101965      | Bottom Rubber Mats for Gen II 96-Well Reactor (101960)  | 25   |
| SHS0034     | Socket Head Screws for Gen II 96-Well and 48-Well Reactors  | 5    |



101968



101967



SHS0034



101480

**Standard 48-Position Photoredox Reaction Block, Gen II for 2mL, 12x32mm Vials**

| Catalog No. | Description  | Qty  |
|-------------|--|------|
| 101480      | 4-Bolt, 48-Well Aluminum Reaction Block for 12x32mm (2mL) Vials. Includes: Each Base Plate, Cover, Screws, PFA Film and Rubber Mat | Each |
| 101968      | Top PFA Films for Gen II 96-Well (101960) and 48-Well (101480) Reactors  | 25   |
| 101967      | Top Rubber Mats for Gen II 96-Well (101960) and 48-Well (101480) Reactors  | 25   |
| SHS0034     | Socket Head Screws for Gen II 96-Well and 48-Well Reactors   | 5    |



101968



101967



SHS0034




101240

**1 Dram 24-Position Photoredox Reaction Block, Gen II with 18mm Well Spacing - for 4mL (1 Dram), 15x45mm Vials**

| Catalog No. | Description   | Qty  |
|-------------|---|------|
| 101240      | 5-Bolt, 24-Well Aluminum Reaction Block for 15x45mm (4mL) Vials. Includes: Base Plate, Cover, Screws, PFA Film and Rubber Mat | Each |
| 101248      | Top PFA Films for Gen II 24-Well Reactor (101240)   | 25   |
| 101247      | Top Rubber Mats for Gen II 24-Well Reactor (101240)   | 25   |
| SHS0114     | Socket Head Screws for Gen II 24-Well Reactors and 96-Well Lightweight Reactors   | 5    |



101248



101247



SHS0114



102960

**Standard 96-Position Photoredox LIGHTWEIGHT Reaction Block, Gen II for 1mL, 8x30mm Glass Inserts**

- 37% lighter than Gen II, 96-well standard 4-bolt reactor
- Great for use on a low-capacity balance with a 500g max limit
- Reduced weight allows for better centrifuge performance
- 9mm Well Spacing

| Catalog No. | Description   | QTY  |
|-------------|---|------|
| 102960      | 4-Bolt, 96-Well LIGHTWEIGHT Aluminum Reaction Block for 1mL, 8x30mm Glass Inserts. Includes Vial Rack, Top Cover, PFA Sheet, Rubber Mats and Screws | Each |
| 102968      | Top PFA Films for Gen II LIGHTWEIGHT 96-Well Reactor (102960)   | 25   |
| 102967      | Top Rubber Mats for Gen II LIGHTWEIGHT 96-Well Reactor (102960)   | 25   |
| 102965      | Bottom Rubber Mat for Gen II LIGHTWEIGHT 96-Well Reactor (102960) - (two required per block)  | Each |
| SHS0114     | Socket Head Screws for Gen II 24-Well Reactors and 96-Well Lightweight Reactors   | 5    |

**96-Well Reactor Assembly Weight\***


Gen II, Lightweight (102960): 477.8g  
- Base only, no lid or screws: 267.8g

Gen II, 4-Bolt (101960): 769.0g  
Gen I, Original (96973): 536.0g


\* Weight without vials. Adding vials will add ~ 131g.



102968



102967



102965



SHS0114

For easy loading, see Assembled 96-Well Vial Trays on page 39

**Parallel Synthesis vs. Photoredox Catalysis Reaction Blocks**

Reactors for both applications are nearly identical. The difference lies in the base (vial rack) of the block. Reactors for parallel synthesis have wells with closed bottoms, whereas photoredox blocks have wells with open holes in the bottom allowing for light transmission via a Lumidox® II Array, as well as to read 2D barcodes.

Parallel Synthesis Block, Gen II



Photoredox Catalysis Block, Gen II



**NEW! Para-dox® Gen II Reaction Blocks for Parallel Synthesis are now available. See website.**



- Validated to have less than 5% solvent loss with prolonged heating above boiling point
- Temperature range -78°C – 150°C
- Quality Stirring with Tumble Stirring
- Silicone Rubber Mats provide compression sealing
- Teflon® PFA Films keep the glass reaction vials from sticking to the silicone rubber mats during heating

### Photoredox Catalysis Reaction Blocks for Screening Chemistry




#### Micro 24-Position Photoredox Reaction Block for 50µL, 4x21mm Micro Vial Glass Inserts

- 4x21mm flat-bottom glass micro vials, 50µL volume (25µL maximum volume if using tumble stir bars)

| Cat. No.    | Description   | Qty  |
|-------------|---|------|
| 24250       | Photoredox 24-Well Micro Block Assembly Includes Vial Rack, Bottom and Top Covers, PFA Sheets, Rubber Mats and Screws | Each |
| 24251       | 24-Well Vial Rack for Photoredox Micro  | Each |
| 24252       | Bottom Plate for Photoredox Micro   | Each |
| 24257       | Bottom PFA Films  | 25   |
| 24258       | Bottom Rubber Mats  | 25   |
| 24256       | Top Cover for 24-Well Photoredox Micro  | Each |
| 24261       | Top PFA Films for 24-Well Photoredox Micro  | 25   |
| 24262       | Top Rubber Mats for 24-Well Photoredox Micro  | 25   |
| 1/2"SCREW96 | 1/2" Bottom Screws  | 100  |
| 1"SCREW     | 1" Top Screws   | 100  |

#### Accessories

| Cat. No.   | Description  | Qty  |
|------------|--|------|
| 10421-Case | Micro Vials - 50µL Flat Bottom Glass Inserts, 4x21mm | 1000 |
| 13257      | Magnetic Tumble Stir Bars, 1.32mm x 1.57mm (Fleas)   | 1000 |
| 24245      | Rack Adapter   | Each |



#### Micro 96-Position Photoredox Reaction Block for 50µL, 4x21mm Micro Vial Glass Inserts

- 4x21mm flat-bottom glass micro vials, 50µL volume (25µL maximum volume if using tumble stir bars)

| Cat. No.    | Description   | Qty  |
|-------------|---|------|
| 96970       | Photoredox 96-Well Micro Block Assembly Includes Vial Rack, Bottom and Top Covers, PFA Sheets, Rubber Mats and Screws | Each |
| 96971       | 96-Well Vial Rack for Photoredox Micro  | Each |
| 96972       | Bottom Plate for Photoredox Micro   | Each |
| 96977       | Bottom PFA Films  | 25   |
| 96978       | Bottom Rubber Mats  | 25   |
| 96976       | Top Cover for 96-Well Photoredox Micro, 6.35mm Thick, 2.54mm Holes  | Each |
| 96981       | Top PFA Films for 96-Well Photoredox Micro  | 25   |
| 96982       | Top Rubber Mats for 96-Well Photoredox Micro  | 25   |
| 1/2"SCREW96 | 1/2" Bottom Screws  | 100  |
| 1"SCREW     | 1" Top Screws   | 100  |

**TECH NOTE:**  
Replace PFA Films with every use for best performance



#### Accessories

| Cat. No. | Description  | Qty  |
|----------|--|------|
| 13257    | Magnetic Tumble Stir Bars, 1.32mm x 1.57mm (Fleas)   | 1000 |
| SD1000   | Milwaukee M4 1/4in. Hex Screwdriver Kit. Includes 2 Batteries and Charger (Note: batteries not available for international shipping) | Each |

- Recommended for top performance!**  
**Assembled 4x21mm Vials in Stackable Tray**
- For easy loading - **Saves Time!**
  - Our thorough QC process ensures **Less Evaporation** compared to loose vials




| Cat. No. | Description   | Qty  |
|----------|---|------|
| 96342    | Assembled Stackable Tray Loaded with 4x21mm Flat-Bottom Vials (10421-CASE), Includes Tray and Vials | Each |





- Validated to have less than 5% solvent loss with prolonged heating above boiling point
- Temperature range -78°C – 150°C
- Quality Stirring with Tumble Stirring
- Silicone Rubber Mats provide compression sealing
- Teflon® PFA Films keep the glass reaction vials from sticking to the silicone rubber mats during heating




### Standard 24-Position Photoredox Reaction Block for 8x30mm, 1mL Glass Inserts



| Cat. No. | Description   | Qty  |
|----------|---|------|
| 24253    | Photoredox 24-Well Block Assembly - Includes Vial Rack, Top Cover, Bottom Cover, PFA Films, Rubber Mats, Screws | Each |
| 24247    | Replacement Top Films for 24 Position, Photoredox Block   | 25   |
| 24248    | Replacement Top Mats for 24 Position, Photoredox Block  | 25   |
| 24259    | Replacement Bottom Films for 24 Position, Photoredox Block  | 25   |
| 24260    | Replacement Bottom Mats for 24 Position, Photoredox Block   | 25   |
| 3/8\"/>  |   |      |

#### Accessories

| Cat. No.   | Description  | Qty  |
|------------|--|------|
| 84001-Case | 1mL Clear Glass Shell Vials, 8 x 30mm  | 1000 |
| 13258      | Stainless Steel Cylinder Stir Bars, 1.98mm x 4.80mm (for 8x30mm Vials)   | 1000 |
| SD1000     | Milwaukee M4 1/4in. Hex Screwdriver Kit. Includes 2 Batteries and Charger (Note: batteries not available for international shipping) | Each |




### Standard 96-Position Photoredox Reaction Block for 8x30mm, 1mL Glass Inserts

| Cat. No. | Description   | Qty  |
|----------|---|------|
| 96973    | Photoredox 96-Well Block Assembly Includes: Vial Rack, Covers, Mats, Films and Screws | Each |
| 96967    | Replacement Top Films for 96-Well Photoredox Block                                    | 25   |
| 96965    | Replacement Top Mats for 96-Well Photoredox Block                                     | 25   |
| 96979    | Replacement Bottom Films for 96-Well Photoredox Block                                 | 25   |
| 96980    | Replacement Bottom Mats for 96-Well Photoredox Block                                  | 25   |
| 1/2\"/>  |   |      |

#### Accessories

#### Recommended for top performance! Assembled 8x30mm Vials in Stackable Tray


- For easy loading - **Saves Time!**
- Our thorough QC process ensures **Less Evaporation** compared to loose vials



| Cat. No. | Description  | Qty  |
|----------|--|------|
| 884001   | Stackable Tray Loaded with 96 8x30mm Shell Vials (84001-CASE), Includes Tray and Vials     | Each |
| 884008   | Stackable Tray Loaded with 96 8x30mm Crimp Top Vials (84008-CASE), Includes Tray and Vials | Each |

See all vial loaders on page 39

| Cat. No. | Description  | Qty  |
|----------|--|------|
| 13258    | SS Cylinder Stir Bars, 1.98mm x 4.80mm (for 8x30mm Vials)  | 1000 |
| SD1000   | Milwaukee M4 1/4in. Hex Screwdriver Kit, includes 2 Batteries and Charger (Note: batteries not available for international shipping) | Each |



#### Top Cover with Wider Holes for 96-Well Reaction Blocks

- Enlarged (4.3mm) Holes for Automated Sampling Instruments
- 316 Stainless Steel

| Cat. No. | Description   | Qty  |
|----------|---|------|
| 96985    | Reaction Block Cover with Enlarged Holes, 6mm Thickness | Each |
| 96986    | Reaction Block Cover with Enlarged Holes, 3mm Thickness | Each |

24245 Rack Adapter (Optional)






13258

#### PRODUCT NOTE:


HTE (High Throughput Equipment) Kits are available (see website)

### HPLC 48-Position Photoredox Reaction Block for 12x32mm, 2mL Glass Vials



Use with 12 x 32mm vials, sold separately.

| Cat. No.    | Description  | Qty  |
|-------------|--|------|
| 48612       | 48-Well <b>Open Top/Open Bottom</b> Aluminum Reaction Block for 12mm (2mL) OD Vials. Includes: Base Plate, Cover, Screws, PFA Film and 2 Rubber Mats | Each |
| 48482       | Rubber Mats for 48 Well (12mm OD) Reaction Blocks  | 25   |
| 48483       | PFA films for 48 Well (12mm OD) Reaction Blocks  | 25   |
| VSCREW48    | 1 1/4" Screws for 48-Well Aluminum Reaction Plate  | 100  |
| 11211-Case  | 11mm, 12 x 32, 2mL Wide Mouth Glass <b>Crimp/Snap</b> Vials  | 1000 |
| 11221-Case  | 11mm, 12 x 32, 2mL Wide Mouth Glass <b>Crimp</b> Vials   | 1000 |
| 11221M-Case | 11mm, 12 x 32, 2mL Wide Mouth Glass <b>Crimp</b> Vials with Marking Spot   | 1000 |




### SS Cylinder Stir Bars for 2mL HPLC Vials

- Parylene Encapsulated
- For Photoredox and Parallel Synthesis Block Systems

| Cat. No. | Description   | Qty  |
|----------|---|------|
| 13258    | Stainless Steel Cylinder Stir Bars, 1.98mm x 4.80mm, for 2mL HPLC Vials | 1000 |

### Magnetic Stir Bars

| Cat. No. | Description                                | Qty |
|----------|--|-----|
| 50225    | Magnetic Stir Bars, 5mm x 2mm, PTFE Coated | 25  |
| 502100   | Magnetic Stir Bars, 5mm x 2mm, PTFE Coated | 100 |



### 48-Position Numbered Crimp Vial Assembly for Reaction Blocks

| Cat. No. | Description   | Qty  |
|----------|---|------|
| 48221M   | 48-Well <b>Numbered</b> Crimp Vial Assembly. Includes: Vials (11221M) Marked 1-48; Stir Bar (13258) in Each Vial; 48-Well Cap Mat (99948) | Each |



### Cap Mat for HPLC (2mL) Vials in 48-Position Block

| Cat. No. | Description  | Qty |
|----------|--|-----|
| 99948    | Clear Silicone/PTFE Cap Mat for HPLC Vials in 48-Position Reaction Block | 5   |



**NEW!**

Gen II Para-dox Reaction blocks are now available. See page 6 for more information.

### 1 Dram 24-Position Photoredox Reaction Block for 15x45mm, 4mL Glass Vials



31531

Use with 15 x 45mm vials, sold separately. See part number 31531-Case.




31554

| Cat. No. | Description   | Qty  |
|----------|---|------|
| 24615    | 24-Well, <b>Open Top/Open Bottom</b> Aluminum Reaction Block for 15mm OD (1 Dram) Vials. 20mm spacing between well centers. Includes: Base Plate, Cover, Screws, PFA Film and 2 Rubber Mats | Each |
| 24282    | High-Temp Silicone/Rubber Cover for 24-Well Reaction Blocks for Vials   | 25   |
| 24283    | PFA Film (0.005" Thick) for 24-Well Reaction Blocks for Vials   | 25   |
| VSCREW24 | 1 3/4" Screws for 24-Well Aluminum Reaction Plate   | 100  |

| Cat. No.   | Description   | Qty  |
|------------|---|------|
| 31531-Case | Advantage™ 13mm, 15 x 45, 4mL Clear Glass Screw Vials, Case | 1000 |
| 31554-Case | PTFE Lined Cap  | 1000 |

### 2 Dram 24-Position Photoredox Reaction Block for 17x60mm, 8mL Glass Vials



31760

Use with 17 x 60mm vials, sold separately. See part number 31760-Case.



31542

| Cat. No. | Description   | Qty  |
|----------|---|------|
| 24617    | 24-Well, <b>Open Top/Open Bottom</b> Aluminum Reaction Block for 17mm OD (2 Dram) Vials. 20mm spacing between well centers. Includes: Base Plate, Cover, Screws, PFA Film and 2 Rubber Mats | Each |
| 24282    | High-Temp Silicone/Rubber Cover for 24-Well Reaction Blocks for Vials   | 25   |
| 24283    | PFA Film (0.005" Thick) for 24-Well Reaction Blocks for Vials   | 25   |
| VSCREW24 | 1 3/4" Screws for 24-Well Aluminum Reaction Plate   | 100  |

| Cat. No.   | Description  | Qty  |
|------------|--|------|
| 31760-Case | Advantage™ 17 x 60mm, 8mL Clear Glass Screw Vials, Case                  | 1000 |
| 31542-CASE | 15mm Solid Black Polypropylene Screw Caps with PTFE/F217 Liners          | 1000 |
| 31543-CASE | Black Open Top PP Cap with PTFE/Silicone (0.065") Septa, 15-425mm Thread | 1000 |





### 24-Position Aluminum Reaction Blocks

- Useful for generating compound arrays (medicinal chemistry)
- Useful for conducting screening reaction conditions
- SBS plate format allows for use in multiple automation applications
- Can be used on tumble stirrers / hot plate stirrers / robotic platforms
- Can easily be used in a glovebox
- Sealed with/ PFA Film and Rubber Mats

### 1 Dram 24-Position Photoredox Reaction Block with 18mm Well Spacing for 15x45mm, 4mL Vials



Use with 15 x 45mm vials, sold separately. See part number 31531-Case.

| Cat. No. | Description   | Qty  |
|----------|---|------|
| 24626    | 24-Well Aluminum Reaction Block with 18mm spacing between well centers. For 15mm OD (1 Dram) Vials. Includes: Base Plate, Cover, Screws, PFA Film and 2 Rubber Mats | Each |
| 24120    | Silicone/Rubber Mats for 24 Well, 18mm Spacing Reaction Blocks  | 25   |
| 24121    | PFA Films for 24 Well, 18mm Spacing Reaction Blocks   | 25   |
| VSCREW24 | 1 3/4" Screws for 24-Well Aluminum Reaction Plate   | 100  |

| Cat. No.   | Description   | Qty  |
|------------|---|------|
| 31531-Case | Advantage™ 13mm, 15 x 45, 4mL Clear Glass Screw Vials, Case | 1000 |
| 31554-Case | PTFE Lined Cap  | 1000 |

### 24-Well Reaction Block with Temperature Transfer Cover

The 24-well, 18mm well spaced Photoredox Reaction Block is also available with a **Temperature Transfer Cover** included. The cover replaces the rubber mat and PFA film, which is not included in this assembly. The cover is also available separately.

| Cat. No. | Description  | Qty  |
|----------|--|------|
| 24627    | 24-Well Reaction Block with 18mm spacing between well centers. For 15mm OD (1 Dram) Vials. Includes: Base Plate, Temperature Transfer Cover and Screws | Each |
| 24122    | Temperature Transfer Cover for 24-Well Reactors with 18mm Well Spacing   | Each |


**NEW!**

### 1 Dram 24-Position Photoredox Reaction Block with 18mm Well Spacing, Gen II for 15x45mm, 4mL Vials



#### Benefits of Gen II Reactors

- Less bolts = less time to assemble
- Less Maintenance
- Port for thermocouple, sensors, etc.
- Corner holes on both the lid and the bottom for accessory attachment
- Longer lasting screws (larger diameter, larger threads)
- Larger holes in lid to accommodate a variety of needle sizes



| Catalog No. | Description   | Qty  |
|-------------|---|------|
| 101240      | 5-Bolt, 24-Well Aluminum Reaction Block for 15x45mm (4mL) Vials. Includes: Base Plate, Cover, Screws, PFA Film and Rubber Mat | Each |
| 101248      | Top PFA Films for Gen II 24-Well Reactor (101240)   | 25   |
| 101247      | Top Rubber Mats for Gen II 24-Well Reactor (101240)   | 25   |
| SHS0114     | Socket Head Screws for Gen II 24-Well Reactors and 96-Well Lightweight Reactors   | 5    |



### Lumidox® II 24-position LED Arrays with 18mm Spacing

Built specifically for use with our **Para-dox® 24-well, 18mm spaced Reactor Blocks**

See page 20 for more information about Lumidox® II


| Wavelength (nm) | Description                                   | Active Cooling Base | Flow-Through Base* | Solid Base*   |
|-----------------|---|---------------------|--------------------|---------------|
|                 |   | Catalog No.         | Catalog No.        | Catalog No.   |
| UV365           | 24-Well LED Array with 18mm Spacing, Lens Mat | LUM22418LA365       | LUM22418LF365      | LUM22418LS365 |
| UV375           | 24-Well LED Array with 18mm Spacing, Lens Mat | LUM22418LA375       | LUM22418LF375      | LUM22418LS375 |
| UV385           | 24-Well LED Array with 18mm Spacing, Lens Mat | LUM22418LA385       | LUM22418LF385      | LUM22418LS385 |
| UV395           | 24-Well LED Array with 18mm Spacing, Lens Mat | LUM22418LA395       | LUM22418LF395      | LUM22418LS395 |
| UV405           | 24-Well LED Array with 18mm Spacing, Lens Mat | LUM22418LA405       | LUM22418LF405      | LUM22418LS405 |
| 420-VIOLET      | 24-Well LED Array with 18mm Spacing, Lens Mat | LUM22418LA420       | LUM22418LF420      | LUM22418LS420 |
| 445-INDIGO      | 24-Well LED Array with 18mm Spacing, Lens Mat | LUM22418LA445       | LUM22418LF445      | LUM22418LS445 |
| 470-BLUE        | 24-Well LED Array with 18mm Spacing, Lens Mat | LUM22418LA470       | LUM22418LF470      | LUM22418LS470 |
| 505-CYAN        | 24-Well LED Array with 18mm Spacing, Lens Mat | LUM22418LA505       | LUM22418LF505      | LUM22418LS505 |
| 527-GREEN       | 24-Well LED Array with 18mm Spacing, Lens Mat | LUM22418LA527       | LUM22418LF527      | LUM22418LS527 |
| 590-AMBER       | 24-Well LED Array with 18mm Spacing, Lens Mat | LUM22418LA590       | LUM22418LF590      | LUM22418LS590 |
| 630-RED         | 24-Well LED Array with 18mm Spacing, Lens Mat | LUM22418LA630       | LUM22418LF630      | LUM22418LS630 |
| WHITE           | 24-Well LED Array with 18mm Spacing, Lens Mat | LUM22418LAWHT       | LUM22418LFWHT      | LUM22418LSWHT |

\* requires the use of a sufficient chiller/cooling device



### Temperature Controlled Reactors (TCR) for 8x30mm (1mL) and 15x45mm (1 Dram) Vials

- Provides extremely uniform thermal control for high throughput experimentation (HTE) with a temperature difference of +/- 1°C
- Capable of enhanced reproducibility in heating and cooling applications
- High quality leak-proof CPC fittings for quick and easy tubing connection and removal
- Compatible with a wide range of heat transfer fluids such as water (down to 5°C), ethylene glycol, polypropylene glycol and silicone-based fluids (ie SYLTHERM™)
- Designed to standard SLAS dimensions (127.75mm x 85.5mm)
- Compatible with auto-samplers and high throughput chemistry devices
- 4mm holes in lid allow for use with most common autosampler needles
- Threaded holes (6-32) in lid for accessory attachment



### Screen & Scale-up with the Same Precision!

#### 24-Position Temperature Controlled Reactor

- For 15x45mm (1 Dram) Vials
- 24 Wells (4 Rows of 6)
- 18mm Well Spacing


| Cat. No.   | Description  | Qty  |
|------------|--|------|
| 249900     | 24-position TCR. Includes: PFA Film, Silicone Rubber Sealing Mat, Rubber Gasket, CPC Fittings, PUR Tubing, Screws. Vials not included. | Each |
| 101247     | Top Rubber Mats for Gen II 24-Position Reactors (with 18mm Spacing)  | 25   |
| 101248     | Top PFA Films for Gen II 24-Position Reactors (with 18mm Spacing)  | 25   |
| SHS0114    | Socket Head Screws for Gen II 24-Position Reactors, 5/16"-18 Threads, 1-1/4" Long  | 5    |
| 31531-Case | 13mm, 15 x 45mm, 4mL (1 Dram) Clear Glass Screw Vials  | 1000 |
| 31554-Case | 13mm Black PP Solid Screw Caps with PTFE/F217 Liners   | 1000 |

#### Lumidox® II 24-position LED Array with 18mm Spacing, Solid Base


- Lens Mat surface
- Requires Thermal Transfer Deck (TTD) connected to chiller

| Cat. No.      | Wavelength (nm) | Cat. No.      | Wavelength (nm) |
|---------------|-----------------|---------------|-----------------|
| LUM22418LS365 | UV365           | LUM22418LS470 | 470-BLUE        |
| LUM22418LS375 | UV375           | LUM22418LS505 | 505-CYAN        |
| LUM22418LS385 | UV385           | LUM22418LS527 | 527-GREEN       |
| LUM22418LS395 | UV395           | LUM22418LS590 | 590-AMBER       |
| LUM22418LS405 | UV405           | LUM22418LS630 | 630-RED         |
| LUM22418LS420 | 420-VIOLET      | LUM22418LSWHT | WHITE           |
| LUM22418LS445 | 445-INDIGO      |               |                 |

#### First generation Photoredox Reactor (uncooled)



#### 48-Position TCR



Simulated heat maps of a standard 96-position Photoredox Reactor Block vs. a 48-position Temperature Controlled Reactor (TCR) when used with a Lumidox® II 48 Position LED Array at full power (stage 5).

#### 24-Position TCR with Temperature Transfer Cover



Perfect for those who prefer to cap their vials instead of sealing them with the standard aluminum lid, rubber mat and PFA film that come with the block assembly. When used in conjunction with a Thermal Transfer Deck+ (TTD+) and external liquid chiller, the cover helps control excessive heat generated by a Lumidox® II LED array and also greatly reduces light bleed.

The Temperature Transfer Cover is also available separately (249913) and can be used with other Para-dox® Gen II 24-position reactor blocks (18mm spacing) for photoredox and parallel synthesis.

| Cat. No. | Description   | Qty  |
|----------|---|------|
| 249800   | 24-Position TCR with Temperature Transfer Cover                           | Each |
| 249913   | Temperature Transfer Cover for Gen II, 24-Position, 18mm Spacing Reactors | Each |

#### Additional Required Parts

| Cat. No.  | Description   | Qty  |
|-----------|---|------|
| 266550    | Para-dox® "Straight-Thru" Thermal Transfer Deck (TTD) | Each |
| 266703    | Daisy Chain Connector for TCR and "Straight-Thru" TTD | Each |
| TCubeEdge | Recirculating Chiller. Operating Range: 0°C - 65°C    | Each |



#### 48-Position Temperature Controlled Reactor

- For 8x30mm (1mL) Vials • 48 Wells (6 Rows of 8) • 9mm Well Spacing

| Cat. No.   | Description   | Qty  |
|------------|---|------|
| 489900     | 48-position TCR. Includes: PFA Film, Blue Fluorosilicone Sealing Mat, Rubber Gasket, Silicone O-Rings, CPC Fittings, PUR Tubing and Screws. | Each |
| 489906     | 1/8" Thick, Blue Fluorosilicone Rubber Sealing Mats for TCR   | 5    |
| 489907     | 0.005" Thick PFA Sealing Films for TCR  | 25   |
| 489908     | 9mm OD, 6mm ID, 1.5mm Wide Silicone O-rings for TCR   | 50   |
| SHS0034    | SS Low-Profile Socket Head Screw w/ Hex Drive, 5/16"-18 Threads, 3/4" Long  | 5    |
| 84001-Case | 1mL Clear Glass Shell Vials, 8 x 30mm   | 1000 |
| 488401     | Vial Tray for TCR, Pre-loaded with 48 8x30 Shell Vials (84001-CASE)   | Ea   |

#### Solid Base Lumidox® II 48-position LED Array for TCR

- Lens Mat surface
- Requires Thermal Transfer Deck (TTD) connected to chiller

| Cat. No.    | Wavelength | Cat. No.    | Wavelength | Cat. No.    | Wavelength |
|-------------|------------|-------------|------------|-------------|------------|
| LUM248LS365 | UV365      | LUM248LS405 | UV405      | LUM248LS527 | 527-GREEN  |
| LUM248LS375 | UV375      | LUM248LS445 | 445-INDIGO | LUM248LS590 | 590-AMBER  |
| LUM248LS385 | UV385      | LUM248LS470 | 470-BLUE   | LUM248LS630 | 630-RED    |
| LUM248LS395 | UV395      | LUM248LS505 | 505-CYAN   | LUM248LSWHT | WHITE      |



### Tech Tip

- Temperature Controlled Reactor (TCR)\* - cools vials
- Thermal Transfer Deck (TTD)\* - cools LED Array (solid base)
- LED Array with solid base - needs TTD and chiller for cooling
- \* Connects to External Liquid Chiller, required



## For Photoredox Catalysis Applications

Taking advantage of recent breakthroughs in LED technology, Analytical now offers the next generation Lumidox® II with higher optical/radiometric power than ever before. Lumidox II offers a multitude of unique and complimentary photonic devices, in varying wavelengths, amplitudes and footprints.

The Lumidox II collection includes new LED Arrays with a wider range of options to choose. Also, check out our LumLamp for experimental applications. All of our Lumidox II illumination devices are controlled by the Lumidox II controller.

### Lumidox II Arrays and LumLamps are available in up to 15 wavelengths:

- 365 (UV365)   405 (UV405)   505 (Cyan)   660 (Deep Red)
- 375 (UV375)   420 (Violet)   527 (Green)   730-IR (Infrared)
- 385 (UV385)   445 (Indigo)   590 (Amber)   White
- 395 (UV395)   470 (Blue)   630 (Red)

Analytical has also developed a myriad of complimentary apparatus that works with Lumidox. These range from SLAS cooling blocks to our new Temperature Controlled Reactors (TCR, page 16). We also offer high throughput electrochemical assemblies that can work with Lumidox to expand your experimentation into electro-photochemical reactions (see page 30).

### Lumidox® II Line of products:

- 96-position LED Arrays
- 48-position LED Arrays for TCR (see page 17)
- 24-position Arrays
- LumLamps
- Lumidox II Controller
- Cell Culture Plate Adapters for Lumidox II LED Arrays (see website)


### Key Features:

- Reduces photochemical reaction time
- Up to 5 selectable output STAGES
- Up to 15 wavelengths available
- Higher optical/radiometric power than previous Lumidox line of products
- Patented LED array design
- Experiment with electro-photochemical reactions by pairing with HTe<sup>-</sup>Chem Electrochemistry Assemblies (see page 30)

Lumidox II Controller



Peripheral Devices



LumLamp  
see page 29

LED Arrays  
see pages 20-24

### Lumidox® II Controller

- Controls the output of a Lumidox Array or LumLamp
- Select one of 5 levels of optical power (STAGES). STAGES are displayed in easy to read, easy to calculate whole numbers.
- Onboard illumination elapsed timer
- Automatic countdown shut-off timer
- **USB Enabled API** - for advanced users running their setups remotely

| Catalog No. | Description  |
|-------------|--|
| LUM2CON     | Lumidox II LED Controller. Includes: Controller, Power Supply, Power Cable, Manual |

### PRODUCT NOTE:

Arrays and LumLamps are tuned in-house. Calibration data is stored onboard (not on the controller), allowing the use of different illuminators with one controller. Note: the controller can operate different Lumidox devices, but only 1 device at a time.

### LED Output STAGES

By default, all Lumidox II devices are factory calibrated with 5 discrete linearly stepped output STAGES. STAGE 1 output is the least radiometric power while STAGE 5 output is the most. STAGES are calibrated to the nearest whole number of radiant flux, and displayed in milliwatts (mW).

#### Linear Step (default)

- STAGE 1 - 50mW
- STAGE 2 - 100mW
- STAGE 3 - 150mW
- STAGE 4 - 200mW
- STAGE 5 - 250mW



### Custom LED Tuned STAGES

Analytical offers custom tuned stages to meet specific and unique requirements. These can range from maximum light output to low output to tight resolution.

**NOTE:** STAGE settings are stored on each individual Lumidox® II peripheral device (LED Array or Lumlamp), and are **not tunable via the controller**. Custom STAGE settings can only be calibrated by Analytical Sales and Services prior to shipment. Please contact Analytical before ordering.

#### High Output (for more aggressive reactions)

- STAGE 1 - 100mW
- STAGE 2 - 200mW
- STAGE 3 - 400mW
- STAGE 4 - 600mW
- STAGE 5 - 800mW

#### High Resolution (for greater precision)

- STAGE 1 - 100mW
- STAGE 2 - 110mW
- STAGE 3 - 115mW
- STAGE 4 - 120mW
- STAGE 5 - 125mW

#### Low Output (good for cell culture and PCR)

- STAGE 1 - 20mW
- STAGE 2 - 40mW
- STAGE 3 - 60mW
- STAGE 4 - 80mW
- STAGE 5 - 100mW

### Lumidox® II LED Arrays

U.S. Patent No. 11,458,447

Our patented LED arrays provide wavelength and power specific illumination to samples for photoredox catalysis applications.

- Up to 15 different color wavelengths are available\*
- 5 output STAGES (customizeable) with calibration data stored on board\*\*
- Available in two Surface Mat styles
- Three options for Base configurations

#### Radiant Flux Values per STAGE - Example

(typical values for 96-Position LED Arrays with Active Cooling Base and Lens Mat)

| Wavelength   | STAGE 1       |           | STAGE 2       |           | STAGE 3       |           | STAGE 4       |           | STAGE 5       |           |
|--------------|---------------|-----------|---------------|-----------|---------------|-----------|---------------|-----------|---------------|-----------|
|              | Per well (mW) | Total (W) | Per well (mW) | Total (W) | Per well (mW) | Total (W) | Per well (mW) | Total (W) | Per well (mW) | Total (W) |
| UV365        | 25            | 2.4       | 55            | 5.3       | 80            | 7.7       | 105           | 10.1      | 135           | 13.0      |
| UV375        | 25            | 2.4       | 50            | 4.8       | 75            | 7.2       | 110           | 10.6      | 140           | 13.4      |
| UV385        | 30            | 2.9       | 65            | 6.2       | 100           | 9.6       | 145           | 13.9      | 165           | 15.8      |
| UV395        | 30            | 2.9       | 65            | 6.2       | 105           | 10.1      | 140           | 13.4      | 170           | 16.3      |
| UV405        | 25            | 2.4       | 60            | 5.8       | 90            | 8.6       | 125           | 12        | 160           | 15.4      |
| 420-Violet   | 30            | 2.9       | 60            | 5.8       | 90            | 8.6       | 120           | 11.5      | 150           | 14.4      |
| 445-Indigo   | 60            | 5.8       | 120           | 11.5      | 195           | 18.7      | 240           | 23.0      | 295           | 28.3      |
| 470-Blue     | 45            | 4.3       | 95            | 9.1       | 140           | 13.4      | 180           | 17.3      | 220           | 21.1      |
| 505-Cyan     | 35            | 3.4       | 65            | 6.2       | 90            | 8.6       | 115           | 11.0      | 135           | 13.0      |
| 527-Green    | 25            | 2.4       | 55            | 5.3       | 80            | 7.7       | 100           | 9.6       | 110           | 10.6      |
| 590-Amber    | 75            | 7.2       | 95            | 9.1       | 115           | 11.0      | 150           | 14.4      | 180           | 17.3      |
| 630-Red      | 30            | 2.9       | 55            | 5.3       | 85            | 8.2       | 115           | 11.0      | 145           | 13.9      |
| 660-Deep Red | 40            | 3.8       | 70            | 6.7       | 105           | 10.1      | 140           | 13.4      | 170           | 16.3      |
| IR730        | 35            | 3.4       | 65            | 6.2       | 95            | 9.1       | 125           | 12.0      | 155           | 14.9      |
| White        | 50            | 4.8       | 100           | 9.6       | 150           | 14.4      | 200           | 19.2      | 300           | 28.8      |


Irradiance charts can be found on our website

\* Most LED arrays (excluding Discovery, page 23) are wavelength-specific - only one wavelength per device  
 \*\* Calibration data is stored onboard (not on the controller), allowing the use of different illuminators with one controller. Controller can only operate 1 device at a time. STAGE settings are NOT user-adjustable and must be configured prior to shipping

### About Array Cooling and Base Options

Lumidox® II LED Arrays can generate a considerable amount of heat at any output stage and therefore need to be cooled. Both **Lens Mat** and **Diffuse Mat** style arrays are available with an **Active Cooling** base, **Solid** base or **Flow-Through** base. Each has its own specific way of cooling the LEDs in the array and can affect different application situations. *Note: solid and flow-through base arrays are not self-cooling and require the use of a sufficient external chiller/cooling device.*

**Active Cooling Base**  
(underside shown)



**Solid Base**  
(shown with Thermal Transfer Deck)



**Flow-Through Base**  
(connects directly to chiller)



### Array Base Comparison

#### Active Cooling Base:

- Self-cooling, no external cooling source required
- Fully conforms to SLAS/ANSI standard dimensions
- Ideal for lower output applications like cell culture and PCR work (may require special adapter)
- Not compatible with tumble stirrers or devices that generate powerful magnetic fields

#### Solid Base:


- External cooling source is required (a Thermal Transfer Deck with a recirculating liquid chiller is recommended)
- A Thermal Transfer Deck (TTD) can be attached to the base when used in conjunction with a recirculating chiller. The TTD can be removed if using a different cooling source, such as a cooling bay
- Fully conforms to SLAS/ANSI standards and can be used with cooling baths, plates, or other cooling chamber (if NOT being used with our Thermal Transfer Deck and recirculating liquid chiller)
- Ideal for applications requiring high output - up to nearly 3x more radiometric power output than Active Base arrays\*

#### NEW! Flow-Through Base:

- Direct connection to recirculating liquid chiller (required). Cooled liquid flows through the array base itself (no need for Thermal Transfer Deck)
- Ideal for applications requiring high output - up to nearly 3x more radiometric power output than Active Base arrays\*
- Light Weight – improved usability with orbital shakers due to lower overall mass
- Shorter overall height than active base array – offers improved compatibility with tumble stirrers (less distance between stirrer and sample)



Active Cooling Base (self-cooling)




Solid Base  
Nearly 3x the radiometric power output of Active Cooling Base arrays (external cooling source required)



Flow-Through Base  
Nearly 3x the radiometric power output of Active Cooling Base arrays (connection to circulating liquid chiller required)

\* For comparison, Radiant Flux Value charts for Active Cooling Base and Solid/Flow-through Base Arrays are posted on our website.  
 Note: All LED arrays can operate in an incubator at 37°C, 95% humidity



LED array with Diffuse Mat under a 250mL flask



LED array with Lens Mat under a Para-dox® 96-well reaction block

### Surface Mat Styles

Array top surfaces come in two styles, **Lens Mat** and **Diffuse Mat**.



#### Lens Mat Surface

- Ultra-clear, molded to fit into the holes of a Para-dox® Reaction Block
- Captures nearly all light emitted by the array's LEDs and directs it into the vials in the reaction block
- Chemically inert silicone



#### Diffuse Mat Surface

- Flat surface, making for easy pairing with SLAS footprint apparatus
- Can be used for niche applications such as illuminating cell culture flasks, reservoir plates, large scale containers, etc.
- Chemically inert silicone

96-Position LED Arrays



96-Position LED Arrays with Lens Mat

| Wavelength (nm) | Description                                    | Active Cooling Base<br>Catalog No. | Flow-Through Base*<br>Catalog No. | Solid Base*<br>Catalog No. |
|-----------------|--|------------------------------------|-----------------------------------|----------------------------|
| UV365           | Lumidox II 96-Position LED Array with Lens Mat | LUM296LA365                        | LUM296LF365                       | LUM296LS365                |
| UV375           | Lumidox II 96-Position LED Array with Lens Mat | LUM296LA375                        | LUM296LF375                       | LUM296LS375                |
| UV385           | Lumidox II 96-Position LED Array with Lens Mat | LUM296LA385                        | LUM296LF385                       | LUM296LS385                |
| UV395           | Lumidox II 96-Position LED Array with Lens Mat | LUM296LA395                        | LUM296LF395                       | LUM296LS395                |
| UV405           | Lumidox II 96-Position LED Array with Lens Mat | LUM296LA405                        | LUM296LF405                       | LUM296LS405                |
| 420-VIOLET      | Lumidox II 96-Position LED Array with Lens Mat | LUM296LA420                        | LUM296LF420                       | LUM296LS420                |
| 445-INDIGO      | Lumidox II 96-Position LED Array with Lens Mat | LUM296LA445                        | LUM296LF445                       | LUM296LS445                |
| 470-BLUE        | Lumidox II 96-Position LED Array with Lens Mat | LUM296LA470                        | LUM296LF470                       | LUM296LS470                |
| 505-CYAN        | Lumidox II 96-Position LED Array with Lens Mat | LUM296LA505                        | LUM296LF505                       | LUM296LS505                |
| 527-GREEN       | Lumidox II 96-Position LED Array with Lens Mat | LUM296LA527                        | LUM296LF527                       | LUM296LS527                |
| 590-AMBER       | Lumidox II 96-Position LED Array with Lens Mat | LUM296LA590                        | LUM296LF590                       | LUM296LS590                |
| 630-RED         | Lumidox II 96-Position LED Array with Lens Mat | LUM296LA630                        | LUM296LF630                       | LUM296LS630                |
| 660-DEEP RED    | Lumidox II 96-Position LED Array with Lens Mat | LUM296LA660                        | LUM296LF660                       | LUM296LS660                |
| 730-IR          | Lumidox II 96-Position LED Array with Lens Mat | LUM296LA730                        | LUM296LF730                       | LUM296LS730                |
| WHITE           | Lumidox II 96-Position LED Array with Lens Mat | LUM296LAWHT                        | LUM296LFWHT                       | LUM296LSWHT                |

\* requires the use of a sufficient chiller/cooling device

See pages 20-21 for details on surface mats and bases



96-Position LED Arrays with Diffuse Mat

| Wavelength (nm) | Description                                       | Active Cooling Base<br>Catalog No. | Flow-Through Base*<br>Catalog No. | Solid Base*<br>Catalog No. |
|-----------------|---|------------------------------------|-----------------------------------|----------------------------|
| UV365           | Lumidox II 96-Position LED Array with Diffuse Mat | LUM296DA365                        | LUM296DF365                       | LUM296DS365                |
| UV375           | Lumidox II 96-Position LED Array with Diffuse Mat | LUM296DA375                        | LUM296DF375                       | LUM296DS375                |
| UV385           | Lumidox II 96-Position LED Array with Diffuse Mat | LUM296DA385                        | LUM296DF385                       | LUM296DS385                |
| UV395           | Lumidox II 96-Position LED Array with Diffuse Mat | LUM296DA395                        | LUM296DF395                       | LUM296DS395                |
| UV405           | Lumidox II 96-Position LED Array with Diffuse Mat | LUM296DA405                        | LUM296DF405                       | LUM296DS405                |
| 420-VIOLET      | Lumidox II 96-Position LED Array with Diffuse Mat | LUM296DA420                        | LUM296DF420                       | LUM296DS420                |
| 445-INDIGO      | Lumidox II 96-Position LED Array with Diffuse Mat | LUM296DA445                        | LUM296DF445                       | LUM296DS445                |
| 470-BLUE        | Lumidox II 96-Position LED Array with Diffuse Mat | LUM296DA470                        | LUM296DF470                       | LUM296DS470                |
| 505-CYAN        | Lumidox II 96-Position LED Array with Diffuse Mat | LUM296DA505                        | LUM296DF505                       | LUM296DS505                |
| 527-GREEN       | Lumidox II 96-Position LED Array with Diffuse Mat | LUM296DA527                        | LUM296DF527                       | LUM296DS527                |
| 590-AMBER       | Lumidox II 96-Position LED Array with Diffuse Mat | LUM296DA590                        | LUM296DF590                       | LUM296DS590                |
| 630-RED         | Lumidox II 96-Position LED Array with Diffuse Mat | LUM296DA630                        | LUM296DF630                       | LUM296DS630                |
| 660-DEEP RED    | Lumidox II 96-Position LED Array with Diffuse Mat | LUM296DA660                        | LUM296DF660                       | LUM296DS660                |
| 730-IR          | Lumidox II 96-Position LED Array with Diffuse Mat | LUM296DA730                        | LUM296DF730                       | LUM296DS730                |
| WHITE           | Lumidox II 96-Position LED Array with Diffuse Mat | LUM296DAWHT                        | LUM296DFWHT                       | LUM296DSWHT                |

\* requires the use of a sufficient chiller/cooling device



Active cooling base shown  
Colors are simulated (UV light is invisible)



Active cooling base shown  
Colors are simulated (UV light is invisible)



Active cooling base shown  
Colors are simulated



Solid base shown\*

Discovery/Screening 96-Position LED Arrays

Lumidox® II Discovery LED Arrays for screening offer scientists an economical and efficient way to evaluate most of the wavelengths that the Lumidox II line of products has to offer. Discovery LED arrays come in varying wavelength combinations and are available with an active cooling, solid or flow-through base.

There are currently 15 wavelength options available in the Lumidox II product line. With Discovery LED arrays, 9 of those wavelengths are combined into 3 units (no need to purchase 9 individual units). The three models of Discovery LED arrays come in these wavelength combinations:

- **Discovery 1 (3 UVs):** 375nm, 385nm, 395nm (4 columns of 8 LEDs / wavelength)
- **Discovery 2 (UV405, Indigo, Blue & Amber):** 405nm, 445nm, 470nm, 590nm (3 columns of 8 LEDs / wavelength)
- **Discovery 3 (Cyan, Green):** 505nm, 527nm (6 columns of 8 LEDs / wavelength)

The six remaining wavelengths in the Lumidox II line - UV365, 420-VIOLET, 630-RED, 660-DEEP RED, 730-IR, and White - are not compatible for placement in Discovery LED arrays due to varying electrical conditions. They can, however, be purchased as single wavelength units to complete your LED array collection.

96-Position Discovery LED Arrays with Lens Mat

| Wavelength (nm)                     | Description  | Active Cooling Base<br>Catalog No. | Flow-Through Base*<br>Catalog No. | Solid Base*<br>Catalog No. |
|-------------------------------------|--|------------------------------------|-----------------------------------|----------------------------|
| UV375                               | Discovery 1 LED Array w/ Lens Mat<br>(4 rows of 8 LEDs per wavelength) | LUM296LAG1                         | LUM296LFG1                        | LUM296LSG1                 |
| UV385                               |  |                                    |                                   |                            |
| UV395                               |  |                                    |                                   |                            |
| UV405                               | Discovery 2 LED Array w/ Lens Mat<br>(3 rows of 8 LEDs per wavelength) | LUM296LAG2                         | LUM296LFG2                        | LUM296LSG2                 |
| 445-INDIGO                          |  |                                    |                                   |                            |
| 470-BLUE                            |  |                                    |                                   |                            |
| 590-AMBER                           | Discovery 3 LED Array w/ Lens Mat<br>(6 rows of 8 LEDs per wavelength) | LUM296LAG3                         | LUM296LFG3                        | LUM296LSG3                 |
| 505-CYAN                            |  |                                    |                                   |                            |
| 527-GREEN                           |  |                                    |                                   |                            |
| <i>Single Wavelength LED Arrays</i> |  |                                    |                                   |                            |
| UV365                               | 96-Position LED Array w/ Lens Mat                                      | LUM296LA365                        | LUM296LF365                       | LUM296LS365                |
| 420-VIOLET                          | 96-Position LED Array w/ Lens Mat                                      | LUM296LA420                        | LUM296LF420                       | LUM296LS420                |
| 630-RED                             | 96-Position LED Array w/ Lens Mat                                      | LUM296LA630                        | LUM296LF630                       | LUM296LS630                |
| 660-DEEP RED                        | 96-Position LED Array w/ Lens Mat                                      | LUM296LA660                        | LUM296LF660                       | LUM296LS660                |
| 730-IR                              | 96-Position LED Array w/ Lens Mat                                      | LUM296LA730                        | LUM296LF730                       | LUM296LS730                |
| WHITE                               | 96-Position LED Array w/ Lens Mat                                      | LUM296LAWHT                        | LUM296LFWHT                       | LUM296LSWHT                |

\* requires the use of a sufficient chiller/cooling device

Thermal Transfer Decks

See page 25 for more information

| Cat. No. | Description   |
|----------|---|
| 266530   | SLAS Thermal Transfer Deck (TTD)                      |
| 266540   | TTD+ (Thermal Transfer Deck with High-end Fittings)   |
| 266550   | SLAS Footprint, "Straight-Thru" Thermal Transfer Deck |





Lumidox II 24-Position LED Arrays are available in the these wavelengths:

- 365 (UV365)
- 395 (UV395)
- 470 (Blue)
- 590 (Amber)
- 375 (UV375)
- 405 (UV405)
- 505 (Cyan)
- 630 (Red)
- 385 (UV385)
- 445 (Indigo)
- 527 (Green)
- White

**Lumidox® II 24-position LED Arrays with 9mm Spacing**

- Typically used with a Para-dox® 24-well Reaction Block for 8x30mm vial inserts (see page 10)
- Can also be used with a HTe<sup>-</sup>Chem Electrochemistry Assembly (see page 30)

| Wavelength (nm) | Description                                  | Active Cooling Base<br>Catalog No. | Flow-Through Base*<br>Catalog No. | Solid Base*<br>Catalog No. |
|-----------------|--|------------------------------------|-----------------------------------|----------------------------|
| UV365           | 24-Well LED Array with 9mm Spacing, Lens Mat | LUM224LA365                        | LUM224LF365                       | LUM224LS365                |
| UV375           | 24-Well LED Array with 9mm Spacing, Lens Mat | LUM224LA375                        | LUM224LF375                       | LUM224LS375                |
| UV385           | 24-Well LED Array with 9mm Spacing, Lens Mat | LUM224LA385                        | LUM224LF385                       | LUM224LS385                |
| UV395           | 24-Well LED Array with 9mm Spacing, Lens Mat | LUM224LA395                        | LUM224LF395                       | LUM224LS395                |
| UV405           | 24-Well LED Array with 9mm Spacing, Lens Mat | LUM224LA405                        | LUM224LF405                       | LUM224LS405                |
| 445-INDIGO      | 24-Well LED Array with 9mm Spacing, Lens Mat | LUM224LA445                        | LUM224LF445                       | LUM224LS445                |
| 470-BLUE        | 24-Well LED Array with 9mm Spacing, Lens Mat | LUM224LA470                        | LUM224LF470                       | LUM224LS470                |
| 505-CYAN        | 24-Well LED Array with 9mm Spacing, Lens Mat | LUM224LA505                        | LUM224LF505                       | LUM224LS505                |
| 527-GREEN       | 24-Well LED Array with 9mm Spacing, Lens Mat | LUM224LA527                        | LUM224LF527                       | LUM224LS527                |
| 590-AMBER       | 24-Well LED Array with 9mm Spacing, Lens Mat | LUM224LA590                        | LUM224LF590                       | LUM224LS590                |
| 630-RED         | 24-Well LED Array with 9mm Spacing, Lens Mat | LUM224LA630                        | LUM224LF630                       | LUM224LS630                |
| WHITE           | 24-Well LED Array with 9mm Spacing, Lens Mat | LUM224LAWHT                        | LUM224LFWHT                       | LUM224LSWHT                |

\* requires the use of a sufficient chiller/cooling device

See pages 20-21 for details on surface mats and bases



**Lumidox® II 24-position LED Arrays with 18mm Spacing**

Built specifically for use with our Para-dox® 24-well, 18mm spaced Reactor Blocks

See page 14 for more information

| Wavelength (nm) | Description                                   | Active Cooling Base<br>Catalog No. | Flow-Through Base*<br>Catalog No. | Solid Base*<br>Catalog No. |
|-----------------|---|------------------------------------|-----------------------------------|----------------------------|
| UV365           | 24-Well LED Array with 18mm Spacing, Lens Mat | LUM22418LA365                      | LUM22418LF365                     | LUM22418LS365              |
| UV375           | 24-Well LED Array with 18mm Spacing, Lens Mat | LUM22418LA375                      | LUM22418LF375                     | LUM22418LS375              |
| UV385           | 24-Well LED Array with 18mm Spacing, Lens Mat | LUM22418LA385                      | LUM22418LF385                     | LUM22418LS385              |
| UV395           | 24-Well LED Array with 18mm Spacing, Lens Mat | LUM22418LA395                      | LUM22418LF395                     | LUM22418LS395              |
| UV405           | 24-Well LED Array with 18mm Spacing, Lens Mat | LUM22418LA405                      | LUM22418LF405                     | LUM22418LS405              |
| 420-VIOLET      | 24-Well LED Array with 18mm Spacing, Lens Mat | LUM22418LA420                      | LUM22418LF420                     | LUM22418LS420              |
| 445-INDIGO      | 24-Well LED Array with 18mm Spacing, Lens Mat | LUM22418LA445                      | LUM22418LF445                     | LUM22418LS445              |
| 470-BLUE        | 24-Well LED Array with 18mm Spacing, Lens Mat | LUM22418LA470                      | LUM22418LF470                     | LUM22418LS470              |
| 505-CYAN        | 24-Well LED Array with 18mm Spacing, Lens Mat | LUM22418LA505                      | LUM22418LF505                     | LUM22418LS505              |
| 527-GREEN       | 24-Well LED Array with 18mm Spacing, Lens Mat | LUM22418LA527                      | LUM22418LF527                     | LUM22418LS527              |
| 590-AMBER       | 24-Well LED Array with 18mm Spacing, Lens Mat | LUM22418LA590                      | LUM22418LF590                     | LUM22418LS590              |
| 630-RED         | 24-Well LED Array with 18mm Spacing, Lens Mat | LUM22418LA630                      | LUM22418LF630                     | LUM22418LS630              |
| WHITE           | 24-Well LED Array with 18mm Spacing, Lens Mat | LUM22418LAWHT                      | LUM22418LFWHT                     | LUM22418LSWHT              |

\* requires the use of a sufficient chiller/cooling device



**TECH NOTE:**

If using in conjunction with a Flow Reactor (page 28), a Thermal Transfer Deck with High-end Fittings (TTD+) must be used



| Cat. No. | Description   | Qty  |
|----------|---|------|
| 266530   | SLAS Thermal Transfer Deck (TTD)  | Each |
| 266540   | TTD+ (Thermal Transfer Deck with High-end Fittings)   | Each |
| 266550   | SLAS Footprint, “Straight-Thru” Thermal Transfer Deck+ (Thermal Transfer Deck with High-end Fittings) | Each |

**Thermal Transfer Decks**

- Must be used with a recirculating liquid chiller or heater
- Can be attached to a Lumidox® II solid base LED arrays for direct cooling of the array


When connected to an external recirculating liquid chiller or heater, our Thermal Transfer Decks can be used as both chilling plates or hot plates. Apply heat or cold directly to your samples or sample vessels (such as Para-dox® Reaction Blocks), or use it to cool your Lumidox® II solid base LED array.

Thermal Transfer Decks may be linked in novel ways to form gradients, or to serve other more experimental needs. Push-to-connect fittings allow for quick tubing insertion/removal. Input and Output ports are user configurable, and may be relocated between side and top faces\*. Units are shipped with side exit ports as standard. Decks are stackable, and multiple decks can be chained together for a total temperature control solution. Compatible with a wide variety of coolants (water, propylene glycol, ethylene glycol, etc.).


\* Ports are not user configurable on the Straight-Thru TTD



Attach a Thermal Transfer Deck+ to a Lumidox II SOLID BASE LED Array for direct cooling of the array



Use an additional TTD+ to transfer heat or cold directly to your reaction block samples



Lumidox / PCR Setup

### Lumidox® II for Cell Culture / Biological Applications

We can support your photon irradiation needs for your cell culture or biological applications, such as:

- Photoconverting live mammalian cells with Dendra protein for fluorescent imaging days after photoconversion
- Testing of light-induced toxicity of certain compounds on different cell lines for photodynamic therapy applications
- Studying the photothermal effects of compounds in live cells
- Testing photodynamic therapy of various innovative compounds on live cells
- The use of photodynamic therapy as a type of cancer treatment

### Lumidox II LED Arrays with Active Cooling Base

Lumidox II LED arrays generate wavelength specific light for your cell culture or biological experimentation.

- Available in 15 wavelengths from UV, Visible through IR.
- Run at 37°C indefinitely (70°C max) at humidity <=99%.
- Active Cooling Base LED arrays are fully compatible with incubation chambers.
- Use **Lens Mat** Arrays for 24 or 96 well flat bottom collection plates, or with Analytical's 48-well TCR.
- Use **Diffused Mat** Arrays for cell culture flasks, petri dishes, larger 12-well and 6-well cell culture plates, etc..
- Active Cooling Base arrays are recommended for bio/cell culture work, but other base configurations are available.

### Cell Culture Plate Adapters


These adapters are designed for **Lens Mat** arrays (*not compatible with diffuse mat arrays*). Plate specific adapters are currently available for 24 and 96-well cell culture plates. They provide consistent and proper alignment of the plate to the LED array. We currently produce adapters for the following plates:

- All our Low Profile collection plates
- Cellvis P96-1.5H-N
- Greiner 655180
- Corning Costar cell culture plates
- Eppendorf cell culture plates (contact us with your specific plate)
- Thermo Scientific Nunc MicroWell plates
- Perkin Elmer ViewPlate (formerly Packard ViewPlate)

Adapters are available for any standard cell culture plate by request.

### Lumidox® II Controller


- Five stages of light intensity to choose from that can be calibrated to your custom levels at the time of purchase - [see STAGE settings, page 19.](#)
- Controller has an onboard timer to reproducibly limit exposure.



Lumidox / Cell Culture Setup



Lumidox II controller (lum2con)



24 well setup with Lumidox 24-well LED Array - See page 24 for array options

### Lumidox® II LED Arrays with Active Cooling Base

Recommended wavelengths for biological and cell culture work. Other wavelengths are available.

| Catalog No. | Description   | Wavelength       |
|-------------|---|------------------|
| LUM296LA445 | Lumidox® II 96-position LED Array with Lens Mat & Active Cooling Base | 445nm - Indigo   |
| LUM296LA470 | Lumidox® II 96-position LED Array with Lens Mat & Active Cooling Base | 470nm - Blue     |
| LUM296LA527 | Lumidox® II 96-position LED Array with Lens Mat & Active Cooling Base | 527nm - Green    |
| LUM296LA630 | Lumidox® II 96-position LED Array with Lens Mat & Active Cooling Base | 630nm - Red      |
| LUM296LA660 | Lumidox® II 96-position LED Array with Lens Mat & Active Cooling Base | 660nm - Deep Red |
| LUM296LAWHT | Lumidox® II 96-position LED Array with Lens Mat & Active Cooling Base | White            |

### Cell Culture Plate Adapters

| Catalog No. | Description  |
|-------------|--|
| A655180     | Adapter for Lumidox II LED Arrays, for use with Greiner 655180 Plate                       |
| AP9615HN    | Adapter for Lumidox II LED Arrays, for use with Cellvis P96-15H-N Well Plates              |
| A204625     | Adapter for Lumidox II LED Arrays, for use with Cell Culture Plate Part # 96624            |
| A204627     | Adapter for Lumidox II LED Arrays, for use with Cell Culture Plate Part # 96626            |
| ACRCLP01    | Adapter for Lumidox II LED Arrays, for use with Corning Cell Culture Plates                |
| A6005430    | Adapter for Lumidox II LED Arrays, for use with Perkin Elmer ViewPlate (Packard ViewPlate) |
| A204628     | Adapter for Lumidox II LED Arrays, for use with Cell Culture Plate Part # 384628           |
| ACRCLP03    | Adapter for Lumidox II LED Arrays, for use with Corning 2592 Plate                         |
| ACRCLP02    | Adapter for Lumidox II LED Arrays, for use with Corning 3631 Plate                         |
| A655097     | Adapter for Lumidox II LED Arrays, for use with Greiner 655097 Plate                       |
| AT11028     | Adapter for Lumidox II 24-position LED Arrays, for use with 24 Deep Well Plate #24128      |

### Cell Culture Plates


| Catalog No. | Description  | QTY |
|-------------|--|-----|
| 96624       | 96-well Cell Culture/Imaging Microplate (Agilent) - Clear with Flat Bottom, 400µL Per Well                   | 40  |
| 96518       | 96 Well, Flat Bottom Clear PS Microplate, 25-340µL working volume per well                                   | 100 |
| 96626       | 96-well Cell Culture/Imaging Microplate (Agilent) - Black with Clear Flat Bottom, 400µL Per Well             | 32  |
| 961515      | Cellvis 96-well Black Polystyrene/Glass Bottom Plate with high performance #1.5 cover glass                  | 20  |
| 96097       | 96 Well, Black/Clear Flat Bottom Well PS Microplate, 25-340µL working volume per well, high binding, sterile | 40  |
| 384628      | 384-well Cell Culture/Imaging Microplate, 140µL Per Well   | 32  |
| 24128       | 24 Deep Well Plate, V-Bottom   | 50  |

Please see our website or catalog for additional collection plates



### PCR Plate Adapter

| Catalog No. | Description  |
|-------------|--|
| AVWRPCR1    | PCR Plate Adapter for Lumidox®II 96-position LED Lens Mat Arrays |



### Cap Mats and Films for PCR Plates

| Cat. No. | Description  | QTY |
|----------|--|-----|
| 96726    | 6mm Silicone/PTFE Cap Mats for 96-well PCR Collection Plates             | 5   |
| 96104    | Pattern Adhesive Resealable Silicone Sealing Film, Round 96-Well Pattern | 5   |



266100  
Flow Reactor

**Flow Reactor**

The Para-dox® Flow Reactor offers an alternative way to conduct photoreactive research on fluidic samples. Used in conjunction with a Lumidox®II Diffuse Mat Array, wavelength-specific light can be applied to sample fluid circulating through 1/8" OD tubing. The Flow Reactor can be liquid cooled or heated and is highly configurable. Add a Thermal Transfer Deck for a complete flow solution.

Users may configure the Flow Reactor in various ways, including multiple input streams using the machined channels. Tubing can be replaced and installed by the user with a tubing ramp tool.

- SLAS/SBS footprint
- Dripless CPC fittings
- 1/8" OD unions for ease of analyte injection
- User-configurable output/input locations, side or top (factory set to side).
- Exposure volume of 2mL using 1/8" OD, 1/16" ID tubing. The exposure volume may be set to 1mL if two streams are used.
- FEP Tubing Provided
- Cradled design reduces light bleed and helps "seat" and align Lumidox unit
- Compatible with water, glycol based coolants, and silicone based coolants


**Flow Reactor**

| Catalog No. | Description  | Qty  |
|-------------|--------------|------|
| 266100      | Flow Reactor | Each |



**Recommended chiller:**


| Cat. No.   | Description        | Qty  |
|------------|--------------------|------|
| TCube Edge | TCube Edge Chiller | Each |



Single Stream Setup



Multistream Setup



Lumidox®II Diffuse Mat Array placement



Flow Reactor (top) paired with a Lumidox®II Diffuse Mat Array (middle) and Para-dox® Thermal Transfer Deck+ (bottom)



*A lamp system specifically designed for your experimental needs!*

**LumLamp**

The LumLamp system stands apart with multiple wavelengths, tight line-widths, and remote (corded) operation. It has radiant flux levels up to 3W. LumLamps are directly controlled with our Lumidox® II Controller, providing optical stability, user configurable timing and output configuration.

With 4 differing reflector types, 3 styles of lens, and up to 12 wavelengths, LumLamp is available in 104+\* possible configurations.

**Lens Types:**



**Reflector Types:**

- 15° Narrow Spot Beam
- 30° Medium Beam
- 40° Wide Beam
- 80° Extra Wide Beam

**LumLamps are available in the following wavelengths:**

- 365 (UV365)
- 375 (UV375)
- 385 (UV385)
- 395 (UV395)
- 405 (UV405)
- 445 (Indigo)
- 470 (Blue)
- 505 (Cyan)
- 527 (Green)
- 590 (Amber)
- 630 (Red)
- White

\* Contact us for custom wavelengths

**Lumidox® II Controller**

- Controls the output of a Lumidox Array or LumLamp
- 5 individually calibrated levels of optical power (STAGES)
- Optical power levels are provided in easy to read, easy to calculate whole numbers
- Onboard illumination elapsed timer
- Automatic countdown shut off timer
- **USB Enabled API** - for advanced users running their setups remotely

**Lumidox II LED Controller**

| Catalog No. | Description   | Qty  |
|-------------|---|------|
| LUM2CON     | Lumidox II LED Controller.<br>Includes: Controller, Power Supply, Power Cable, Manual | Each |





Lumidox Controllers come with one pair of UV protective goggles. Additional goggles can be purchased.



**UV Safety Goggles**

| Catalog No. | Description                  | Qty  |
|-------------|------------------------------|------|
| UVGoggles   | Additional UV Safety Goggles | Each |





### HTE<sup>-</sup>Chem Electrochemistry Assemblies

Analytical Sales' new HTE<sup>-</sup>Chem line of products allows for broad electrochemical reactions (such as electrosynthesis, organic electrochemistry, electrophotocatalysis, etc.) to be carried out on Analytical Sales standard 24-well minirack platform. HTE<sup>-</sup>Chem allows you to accelerate your electrochemistry workflow by permitting multiple simultaneous constant current, constant voltage, and/or electrophotocatalytic experiments to be run.

There are two optimized setups to choose from:

The **constant current** setup allows the user to set a precise current, which will be maintained by the power supply throughout the reaction. The voltage will automatically be varied based on the changing electrical conditions as the reaction progresses.

The **constant voltage** setup allows the user to set a precise voltage, which will be maintained by the power supply throughout the reaction. In this setup, the current will automatically be varied based on the changing electrical conditions as the reaction progresses.

### HTE<sup>-</sup>Chem Assemblies & Kits




| Cat. No. | Description  | Qty  |
|----------|--|------|
| 700100   | Constant Current Electrochemistry Assembly                     | Each |
| 700150   | Constant Current Electrochemistry Kit (Includes Electrode Kit) | Each |
| 700200   | Constant Voltage Electrochemistry Assembly                     | Each |
| 700250   | Constant Voltage Electrochemistry Kit (Includes Electrode Kit) | Each |

### Electrodes for HTE<sup>-</sup>Chem Assemblies

Analytical Sales offers a myriad of electrode choices which can be used as either cathodes or anodes depending on your experimental needs.


| Cat. No.   | Description   | Qty  |
|------------|---|------|
| 700500     | Graphite Electrodes   | 30   |
| 700500-250 | Graphite Electrodes - Bulk Pack   | 250  |
| 700550     | Tin Electrodes  | 25   |
| 700575     | Cadmium Electrodes  | 25   |
| 700600     | Zinc Electrodes   | 25   |
| 700600-250 | Zinc Electrodes - Bulk Pack   | 250  |
| 700650     | Platinum Electrodes   | 25   |
| 700656     | Platinum Electrodes - Row Pack  | 6    |
| 700675     | Iron Electrodes   | 25   |
| 700700     | Stainless Steel Electrodes  | 25   |
| 700700-250 | Stainless Steel Electrodes - Bulk Pack  | 250  |
| 700725     | Cobalt Electrodes   | 25   |
| 700750     | Nickel Electrodes   | 25   |
| 700800     | NiChrom Electrodes  | 25   |
| 700800-250 | NiChrom Electrodes - Bulk Pack  | 250  |
| 700850     | Copper Electrodes   | 25   |
| 700900     | Aluminum Electrodes   | 25   |
| 700950     | Magnesium Electrodes  | 25   |
| 700400     | Electrode Kit. Includes One Pack of Each of Graphite, Zinc, Stainless Steel, Nickel, NiChrom, Aluminum, Copper and Magnesium Electrodes (platinum not included) | Each |

Check website for the latest selection of electrodes



### Power Supply and Accessories

| Cat. No.   | Description   | Qty  |
|------------|---|------|
| 700050     | Calibrated DC Power Supply, 4 Output Multi Range, 420W Combined Output              | Each |
| 700060     | Breakout Cable Assembly, for use with Calibrated DC Power Supply                    | Each |
| 700070     | Ribbon Cable (if using alternate power supply)                                      | Each |
| 701020     | Replacement PCB Jumpers (for constant current setup)                                | 100  |
| 701015     | Upper Gasket (recommended: replace after 3-5 uses)                                  | 5    |
| 701005     | Lower Gasket (recommended: replace after 3-5 uses)                                  | 5    |
| 701025     | PFA Sheet for HTE <sup>-</sup> Chem Assembly (recommended: replace after every use) | 25   |
| 24260      | Bottom Rubber Mats  | 25   |
| 84001-CASE | 1mL Clear Glass Shell Vial  | 1000 |



For experimentation with electrophotocatalytic reactions, HTE<sup>-</sup>Chem assemblies are fully compatible for use with Lumidox®II 24-well, 9mm spaced LED arrays.



See page 24 for more information

### Flow Electrolysis

Analytical is pleased to offer a solution for electrolysis experimentation in the form of a standardized flow electrolysis setup. This new and upcoming technology allows for rapid and efficient production of complex chemical compounds via flow electrolysis. With an interelectrode gap smaller than a millimeter, and the introduction of flow geometries that promote turbulent mixing, our Flow Electrolysis Cell allows for heightened selectivity and yield when compared to flask based reactors. Additionally, the Flow Electrolysis Cell requires minimal new equipment to use. With fittings for 1/8" OD HPLC tubing, fluid flow and electrodes can be connected to power supplies via standard connectors will automatically be varied based on the changing electrical conditions as the reaction progresses.

#### Undivided Cell (single stream)

In an undivided cell, both electrodes are washed with the same electrolyte and only one fluid circuit is used; the opposite side of the cell is sealed with HPLC plugs and an electrode without through-holes. This allows for low-resistance reactions which will output combined products.



#### Flow Electrolysis Undivided Cell (single stream) Assembly

##### Assembly includes:

- Impervious Graphite Plate Electrode (anode)
- 316L Stainless Steel Plate Electrode (cathode)
- Flow Channel PFA Films (1 Meandering, 1 Tangential, 1 Simple, 1 Fin Separator)
- Isolation Layer PFA Films (2)
- Orange Silicone Rubber Alignment Gaskets (2)
- Viton™ Chemical-resistant Internal O-rings, 2mm ID (2)
- Viton™ Chemical-resistant External O-rings, 3mm ID (2)
- Spring Test Probes/Pogo Pins (1 red, 1 black)
- Aluminum End Plates (2)
- Fittings and Hardware


| Cat. No. | Description   | Qty  |
|----------|---|------|
| F1088100 | Flow Electrolysis Undivided Cell Assembly (single stream) | Each |

##### Channel Film Volume

| Film Type     | Surface Area (mm <sup>2</sup> ) | Volume at 0.254mm Thick (mm <sup>3</sup> ) | Volume at 0.127mm Thick (mm <sup>3</sup> ) |
|---------------|---------------------------------|--|--|
| Simple        | 681.06193                       | 0.26813                                    | 0.13407                                    |
| Tangential    | 485.74568                       | 0.19124                                    | 0.09562                                    |
| Meandering    | 631.43901                       | 0.24860                                    | 0.12430                                    |
| Fin Separator | 474.58642                       | 0.18685                                    | 0.09342                                    |

#### Divided Cell (dual stream)

In a divided cell the electrodes and their flow channels are separated by a Nafion membrane, allowing for reactions to be performed that produce separate output streams. This can be advantageous if the two electrolyte solutions are incompatible or if the products of their reaction are troublesome to separate.





#### Flow Electrolysis Divided Cell (dual stream) Assembly

##### Assembly includes:

- Impervious Graphite Plate Electrode with holes (anode)
- 316L Stainless Steel Plate Electrode (cathode)
- Flow Channel PFA Films (2 Meandering, 2 Tangential, 2 Simple, 2 Fin separator)
- Isolation Layer PFA Films (2)
- Nafion Ion Exchange Channel Separator Membrane
- Orange Silicone Rubber Alignment Gaskets (2)
- Viton Chemical-resistant Internal O-rings, 2mm ID (4)
- Viton Chemical-resistant External O-rings, 3mm ID (4)
- Spring Test Probes/Pogo Pins (1 red, 1 black)
- Aluminum End Plates (2)
- Fittings and Hardware


| Cat. No. | Description   | Qty  |
|----------|---|------|
| F1088200 | Flow Electrolysis Divided Cell Assembly (dual stream) | Each |



#### Electrode Plates and Consumable Accessories

| Cat. No. | Description  | Qty  |
|----------|--|------|
| F1088109 | 316L Stainless Steel Plate Electrode (cathode)       | Each |
| F1088119 | Titanium Electrode (cathode)                         | Each |
| F1088108 | Impervious Graphite Plate Electrode (anode)          | Each |
| F1088208 | Impervious Graphite Plate Electrode (anode) w/ Holes | Each |
| F1088405 | PFA Isolation Layer Films                            | 25   |
| F1088504 | Silicone Rubber Alignment Gaskets                    | 10   |
| F1088406 | PFA Simple Channel Films                             | 25   |
| F1088403 | PFA Meandering Channel Films                         | 25   |
| F1088407 | PFA Tangential Mixer Channel Films                   | 25   |
| F1088410 | PFA Fin Separator Channel Films                      | 25   |
| F1088201 | Nafion Ion Exchange Channel Separator                | Each |
| F1088302 | Viton Chemical-resistant External O-rings            | 50   |
| F1088301 | Viton Chemical-resistant Internal O-rings            | 50   |
| F1088303 | Spring Test Probes / Pogo Pins                       | 10   |
| F1300050 | Flow Electrolysis Power Supply                       | Each |

**EquaVAP**<sup>®</sup>  
Blowdown Evaporators



EquaVAP (23096) shown with Gen II Para-dox<sup>®</sup> Reaction Block (101960)

### EquaVAP<sup>®</sup> 96-Well Blowdown Evaporators - for Reaction Blocks

**Rapidly evaporates common organic solvents**

- **Internal Flow Equalizers** - Distributes equal output across all needle ports to ensure even and symmetrical evaporation.
- **Stainless Steel Needle Tips** - Long term robust use, sturdy, wide solvent compatibilities
- **Step-down Height Adjustment** - Offers positive repeatable height positions
- **Stand Alone** - Simple to operate, connects to Nitrogen or air supply with a standard socket
- **Small Footprint (8" base)** - Fits into small glove boxes & standard fume hoods



### 96-Well EquaVAP Evaporator

- Use with 96-well Reaction Blocks with 1mL Vials or 96-well Polypropylene Collection Plates

| Cat. No. | Description                         |
|----------|-------------------------------------|
| 23096    | EquaVAP 96-Well Blowdown Evaporator |


#### EVAPORATION TIMES

| Chemical                | Avg Evap Time* (HH:MM:SS) |
|-------------------------|---------------------------|
| Dichloromethane 99.9%   | 00:00:23                  |
| Acetone 100%            | 00:02:04                  |
| Methanol 99.9%          | 00:06:09                  |
| Ethyl Acetate           | 00:03:04                  |
| Acetonitrile            | 00:06:37                  |
| Isopropyl Alcohol 99.0% | 00:08:48                  |
| Water (DI)              | 01:23:46                  |

**Conditions: 100 uL of solvent evaporated at AIR flow of 100 L/min (using 1mL round deep well plate)**

\*Individual results may vary

### EquaVAP<sup>®</sup> 48 and 24-Well Blowdown Evaporators - for Reaction Blocks



### 48-Well EquaVAP Evaporator

- Use with 48-well Aluminum Reaction Blocks with 2mL Vials (101480), (48612), (48012)

| Cat. No. | Description  |
|----------|--|
| 23048    | EquaVAP 48-Well Blowdown Evaporator for 48-Well Reaction Blocks with 2mL vials |



### 24-Well EquaVAP for 9mm Vial Spacing

- Use with 24-well optimization blocks with 9mm vial spacing
- 24-well block fits into 96-well footprint with use of rack adapter (24245)

| Cat. No. | Description  |
|----------|--|
| 23024-09 | EquaVAP 24-Well Evaporator for 24-Well Reaction Blocks with 9mm Vial Spacing |



### 24-Well EquaVAP Evaporator (20mm)

- 20mm center-to-center needle distance
- Use with 24-well Aluminum Reaction Blocks with 1 Dram (4mL) or 2 Dram (8mL) Vials (24015), (24017), (24615), (24617)

| Cat. No. | Description  |
|----------|--|
| 23024-20 | EquaVAP 24-Well Evaporator for 24-Well Reaction Blocks with 20mm Spacing |



### 24-Well EquaVAP Evaporator (18mm)

- 18mm center-to-center needle distance
- Use with 24-well Aluminum Reaction Blocks with 18mm Well Spacing (101240), (24626)

| Cat. No. | Description  |
|----------|--|
| 23024-20 | EquaVAP 24-Well Evaporator for 24-Well Reaction Blocks with 18mm Spacing |

Shown with Gen II Para-dox<sup>®</sup> Reaction Block (101240)

#### PRODUCT NOTE:

Requires regulated air supply set between 80-110 psi




**Advantage™ 1mL Collection Plates**  
Round Well with **Round Bottom**

| Cat. No. | Description  | Qty |
|----------|--|-----|
| 17P687Z  | 1mL 96-Well Collection Plate with Round Well Bottoms | 20  |
| 17P687   | 1mL 96-Well Collection Plate with Round Well Bottoms | 20  |




**Advantage™ 2mL 96-Well Collection Plates**  
Square Well with **V-Shaped** or **Round Bottom**

| Cat. No. | Description   | Qty |
|----------|---|-----|
| 27P687   | 2mL 96-Well Collection Plate with <b>Round-Bottom</b> Wells | 20  |
| 59623-23 | 2mL 96-Well Collection Plate with <b>V-Bottom</b> Wells     | 10  |




**Advantage™ 96-Well Cap Mats - Round Well**

| Cat. No. | Description                                  | Qty |
|----------|--|-----|
| 965075   | Purple Pre-Scored Ultra Thin Round Cap Mat   | 5   |
| 964075   | Clear Pre-Slit Ultra Thin Round Well Cap Mat | 5   |
| 96057    | Autosampler Compatible 1mL Cap Mat with "X"  | 20  |



**Advantage™ 96-Well Cap Mats - Square Well**

| Cat. No. | Description                                      | Qty |
|----------|--|-----|
| 965085   | Purple Ultra Thin <b>Pre-Slit</b> Square Cap Mat | 5   |
| 964008   | Mighty Mat Purple Ultra Thin Square Cap Mat      | 5   |
| 2726     | 2mL Pierceable Cap Mat with "X" on Surface       | 20  |



**24-Well Collection Plates and Cap Mats**

| Cat. No. | Description  | Qty |
|----------|--|-----|
| 24727    | 24-Well Collection Plate with Round-Bottom Wells. 10mL/well, 240mL max               | 25  |
| 24108    | 24-Well Collection Plate with Pyramid-Bottom Wells. 10mL/well, 240mL max             | 25  |
| 24777    | Silicone/PTFE <b>Cap Mat</b> for 24-well Collection Plates                           | 5   |
| 24128    | 24 Well Plate, V-Bottom, 15mL Per Well, Polypropylene (127.2L x 85.3W x 63.2H in mm) | 50  |

*Fits nicely with Lumidox® II 24-well LED arrays*





**Thermal Adapter Plates**

- For improved temperature transfer to samples
- For use with matching polypropylene collection plates

Our Thermal Adapter Plates conform to the well shape of their complementing polypropylene collection plate. They ensure uniform thermal transfer for heating and cooling when used in conjunction with a heating block or chiller (not included). The plate features a SLAS footprint and easy to identify purple anodized coating.

*Analytical Sales & Services can make Thermal Adapter Plates to accommodate any shape and style polypropylene collection plate we offer. Please call 973-616-0700 for more information.*



**Thermal Adapter Plates**


| Cat. No.    | Description  | Qty  |
|-------------|--|------|
| 967720TAP   | Thermal Adapter Plate for 967720 PP Collection Plate     | Each |
| 967720      | 2mL Round Well/Round Bottom 96-Well Collection Plate     | 25   |
| 59623-23TAP | Thermal Adapter Plate for 59623-23 PP Collection Plate   | Each |
| 59623-23    | 2mL Round Well/V-Bottom 96-Well Collection Plate         | 25   |
| 968810TAP   | Thermal Adapter Plate for 1mL TrueTaper Collection Plate | Each |
| 968810      | 1mL TrueTaper Collection Plate                           | 25   |
| 968820TAP   | Thermal Adapter Plate for 2mL TrueTaper Collection Plate | Each |
| 968820      | 2mL TrueTaper Collection Plate                           | 25   |
| 17P687ZTAP  | Thermal Adapter Plate for 17P687Z PP Collection Plate    | Each |
| 17P687Z     | 1mL 17P687Z PP Collection Plate                          | 25   |

**Vacuum Manifold Filtration System for SPE Sample Preparation**

- Sturdy, clear acrylic construction
- Adjustable vacuum control
- Easy access to filtrate
- Compatible with all standard filter bottom microplates
- Collect into any storage plate
- Compatible with robotic handling

The Advantage™ Vacuum Manifold Filtration System allows you to collect into any storage plate (including low volume collection plates) and is compatible with robotic handling systems. It uses an integral flat gasket between the collection plate and the filtration plate which stays totally secure within the system. A spacer is supplied for use with low volume collection plates. The Vacuum Manifold is constructed of tough acrylic for optimum performance and visual accessibility.

| Cat. No. | Description   | Qty  |
|----------|---|------|
| 96844    | Vacuum Manifold Filtration System for SPE Sample Filtration | Each |



**Drug Discovery 96-Well Filter Plates**

- Low Protein Binding
- General Filtration

| Cat. No. | Description                  | Qty |
|----------|------------------------------|-----|
| 964PP45  | 400µL Hydrophobic PP, 0.45µm | 25  |

**General Filtration 96-well Filter Plates**

- After SPE
- Dilute and Shoot
- Use prior to LCMS and Microarraying

*Low Protein Binding - Sample in Solvent*

| Cat. No. | Description                | Qty |
|----------|----------------------------|-----|
| 96245-10 | 2mL Hydrophilic PP, 0.45µm | 10  |
| 96254-10 | 2mL Hydrophobic, PP 0.45µm | 10  |

**Diatomaceous 96-Well Filter Plates**


- pH 1-13
- No Pre-Treatment of the Bed is Necessary
- Packed with Flux Calcinated Diatomaceous Earth

| Cat. No. | Description                   | Qty |
|----------|-------------------------------|-----|
| 96160-5  | 1mL Diatomaceous Filter Plate | 5   |
| 96260-5  | 2mL Diatomaceous Filter Plate | 5   |



**General Filtration 384-well Filter Plates**

| Cat. No. | Description                                       | Qty  |
|----------|---|------|
| 38407    | 384-Well Filter Plates, 140µL, Glass Fiber, 0.7µm | 10   |
| 384603   | Poly DVB SPE, 384 Well filter plate, 3mg/well     | Each |



**24-well Filter Plates**

| Cat. No.   | Description   | Qty  |
|------------|---|------|
| 24152      | 24 Well Filter Plates, PE 25µm, 15mL                  | 15   |
| 2415SCX500 | 24 well SCX (Strong Cation exchange) 15 mL/500 mg bed | Each |
| 241010     | 24 well Filter plates, Hydrophobic PP, 10 µm          | 25   |

**Assembled Vials in Stackable Trays**

- For easy loading - **Saves Time!**
- Our thorough QC process ensures **Less Evaporation** compared to loose vials

*Recommended for top performance!*


**8x30mm Crimp Top Vials in Loaders**

**8x30mm Crimp Top Vials for TCR (48)**



| Cat. No. | Description  | Qty  |
|----------|--|------|
| 488408   | Assembled Stackable Tray Loaded with 48 8x30mm Crimp Top Vials (84008-CASE) for TCR, Includes Tray and Vials | Each |

**8x30mm Crimp Top Vials (96)**



| Cat. No. | Description  | Qty  |
|----------|--|------|
| 884008   | Assembled Stackable Tray Loaded with 96 8x30mm Crimp Top Vials (84008-CASE), Includes Tray and Vials | Each |


**Shell Vials in Loaders**

**8x30mm Glass Inserts for TCR (48)**



| Cat. No. | Description  | Qty  |
|----------|--|------|
| 488401   | Assembled Stackable Tray Loaded with 48 8x30mm Flat-Bottom Vials (84001-CASE) for TCR, Includes Tray and Vials | Each |

**8x30mm Glass Inserts (96)**



| Cat. No. | Description   | Qty  |
|----------|---|------|
| 884001   | Assembled Stackable Tray Loaded with 8x30mm Shell Vials (84001-CASE), Includes Tray and Vials | Each |

**4x21mm Vials (96)**



| Cat. No. | Description   | Qty  |
|----------|---|------|
| 96342    | Assembled Stackable Tray Loaded with 4x21mm Flat-Bottom Vials (10421-CASE), Includes Tray and Vials | Each |

**5x31mm Vials (96)**





| Cat. No. | Description   | Qty  |
|----------|---|------|
| 96242    | Assembled Stackable Tray Loaded with 5x31mm Flat-Bottom Vials (20303-CASE), Includes Tray and Vials | Each |


**Powder Transfer Plate for 4x21mm Vials**

*Easily dispense powder into 4x21mm (50µL) Micro Vials*

- Wide funneled holes guide powder into vials
- Grooves on bottom surface "lock" vials into position to avoid spillage



- Wide funneled holes guide powder into vials
- Grooves on bottom surface "lock" vials into position to avoid spillage



1. Lower plate onto vials and let the vials "settle" into the grooves on the bottom surface of transfer plate
2. Powder can now be easily dispensed into each vial


**Use with:**

- 96-Well Polypropylene Tray with 4x21mm Vials (96342)
- 4x21mm Vials (10421-Case) in 96-Well Micro Optimization Block (96970)

| Cat. No. | Description                          | Qty  |
|----------|--------------------------------------|------|
| 96T100   | Powder Transfer Plate for 4x21 Vials | Each |

**Powder Transfer Plate for 8x30mm Vials**

*Use to transfer powder from vial to vial or to mix powders between vials*



Use with 96-Well Stackable Trays with 8x30mm Vials (884001 Shell), (884008 Crimp)

| Cat. No. | Description  | Qty  |
|----------|--|------|
| 96T105   | Powder Transfer Plate for 8x30mm Shell Vials             | Each |
| 884001   | 96-Well Stackable Tray Loaded with 8x30mm Shell Vials    | Each |
| 96T102   | <b>NEW!</b> Powder Transfer Plate for 8x30mm Crimp Vials | Each |
| 884008   | 96-Well Stackable Tray Loaded with 8x30mm Crimp Vials    | Each |

Flip repeatedly to mix powders between vials

**We are here to serve your needs quickly and efficiently! How do we do it?**

- Personalized consultations and direct technical support
- Immediate response times
- Various convenient methods of contact
  - Form submissions through our website
  - Direct email (sales@analytical-sales.com)
  - Call us, a real person will answer!
- Virtual meetings with our sales team
  - Same day availability
  - Flexible scheduling
  - Accommodating international time zones
- Exceptional customer service
  - Receive a quote within less than one day
  - Email responses within hours, often minutes (during normal hours of operation)

**How to Order**

To place your order, you will need to supply the catalog number, a brief description and size, or the particular specifications when indicated.

**You can order by:**

- **Phone:** 973-616-0700
- **Fax:** 973-616-0133
- **Email:** orders@analytical-sales.com
- **Online store:** www.analytical-sales.com

**Terms of Payment**

If you have an account with Analytical Sales, we will bill you for your purchases. All prices are F.O.B. Flanders, NJ. Terms of payment are net 30 days. To open an account, please call us. We also accept Visa®, Mastercard®, and American Express®. PayPal® is also available when ordering from our website.

**Shipping**

All items will be shipped via FedEx® Ground or common carrier unless otherwise instructed. Please examine all items immediately upon receipt. If you notice that an item was damaged in transit, it's important that you get a "damage notation" from the driver. If you notice damage upon unpacking an item, be sure to save all containers and packing materials. Please notify us immediately for instructions.

**Returns**

To receive credit for any product you return, you must first receive authorization. Please contact us for instructions. Returns must be made within 10 days of receiving authorization.

**Pricing**

To see pricing please visit www.analytical-sales.com/catalog/pricing. All prices are subject to change without notice.



179 Rt 206 • Flanders, NJ • 07836

Phone: 973-616-0700 • Fax: 973-616-0133

Email: [info@analytical-sales.com](mailto:info@analytical-sales.com)

Website: [www.analytical-sales.com](http://www.analytical-sales.com)

