


Lumidox II LED Arrays are robust devices that will provide countless hours of operation when used and maintained properly. When using in an environment where volatile organic compounds (VOCs) are present, device longevity and expected operation are dependent upon proper workspace ventilation. Lack of ventilation can cause VOC contamination and possible damage to your device. However, VOC damage is usually temporary, and a procedure for restoring your device is explained at the bottom of this page.

Potential LED Damage due to VOC Contamination

Contaminants can darken LEDs, resulting in rapid light output degradation. This effect can be amplified when exposure to environmental pollutants is combined with high temperatures and/or high intensities of high-energy photons in an environment with little to no ventilation.

When contamination occurs, it darkens the die surface, blocks the light, and thus dramatically depreciates the LED's luminous flux. Some materials may be more likely to cause LED darkening under specific conditions. Contaminants can also damage LEDs through corrosion, causing open circuits. The silicone encapsulant can be damaged due to swelling, bulk discoloration, and loss of adhesion.



VOC Concentration – High concentrations of VOCs will result in more severe darkening. The threshold levels are not well established. System designs that trap VOCs exacerbate the effect.

VOC-induced LED darkening results in discoloration, causing the LED surface to take on a yellow or brown color. Darkening can dramatically depreciate the luminous flux performance of the LED, causing it to emit significantly less light. With white LEDs, the color point can be shifted due to VOC darkening.

Some contaminants that may cause VOC Contamination

Note: The concentration levels that result in degradation are not easily quantified.

Acetic Acid	Cutting Fluids	Iodine	Phenols
Acetone	Cyanoacrylates	Isoamyl Alcohol	Polymer Plasticizers
Acrylates (all types)	Cyclo Hexane	Many types of inks (especially Sharpie® pens)	Potassium Hydroxide
Ammonia	Diamyl Phthalate	Kapton tape with acrylic-based adhesives (silicone-based adhesives are Ok)	Rubber
Ammonium Sulfide	Dichloromethane	Methacrylic Acid-Methyl Ester (MMA)	Silicone Oil
Aniline	Dipropylene Glycol Monomethyl Ether	Methanol	Silicone Rubber (Acetic Acid Curing)
Benzene	Epichlorohydrin	Methyl Amine	Sodium Hydroxide
Bleach	Ethyl Acetate	Methyl Ethyl Ketone (MEK)	Sodium Sulfide
Bromine	Ethyl Amine	Methyl Isobutyl Ketone	Sulfur-Containing Compounds
Butadiene Rubber	Ethyl Bromide	Methyl or Ethyl Acetates	Sulfuric Acid
Butyl Rubber	Ethylene Chlorohydrin	Methylene Chloride	Sulfur Dioxide
Camphor Oil	Formaldehyde	Nitric Acid	Sulphuryl Chloride
Carbolic Acid	Glycol Ethers	Pesticides	Toluene
Carbon Disulfide	Halogenated Hydrocarbons	Petroleum products (all)	Trimethyl Hexamethylene Diamine
Carbon Tetrachloride	Hydrochloric Acid		Xylene
Chlorine-Containing Compounds	Hydrofluoric Acid		
Chlorobutyl Rubber	Hydrogen Peroxide		
Chlorosulphonated Polyethylene			

VOC-induced darkening is often reversible. The silicone dome and LED chip are generally not permanently affected by VOC contamination and contaminants can be driven out of the package by using a burn-in process. **In an open area with adequate ventilation, run your array on stage 5 for 2-3 hours with a fan blowing on the LEDs to completely burn away all the VOCs and revert the array back to its original state.** VOC-driven LED darkening will likely not be severe and will dissipate quickly.