



## 20mm Crimp Style Head Space Vials

### Crimp style vials, made of borosilicate glass

For continuous results it is very important the head space vial and cap form a tight seal. In this type of analysis, the sample is sealed and heated to very high temperatures so the concentration in the gas phase (head space) is equal to the amount in the liquid phase. An aliquot from the headspace is then analyzed.



Cat. No.	Description	Qty
62321-3	6mL Clear Glass Crimp Vial, Bevel Top / Flat Bottom, 22 x 38mm	100
62101-3	10mL Clear Glass Crimp Vial, Bevel Top / Flat Bottom, 23 x 46mm	100
62201-3	20mL Clear Glass Crimp Vial, Bevel Top / Flat Bottom, 23 x 75mm	100

#### Recommended for PE Autosamplers

67220-3	20mL <b>Amber</b> Glass Crimp Vial, Bevel Top / Round Bottom, 23 x 75mm	100
---------	---	-----

#### Recommended for CTC Autosamplers

62282-3	6mL Clear Glass Crimp Vial, Beveled Top / Round Bottom, 22 x 38mm	100
61182-3	10mL <b>Amber</b> Glass Crimp Vial, Beveled Top / Round Bottom, 23 x 46mm	100
61182-3-Case	10mL <b>Amber</b> Glass Crimp Vial, Beveled Top / Round Bottom, 23 x 46mm	1000



## 20mm Crimp Caps with Liners

### Available in Aluminum, Aluminum/Tin and Magnetic Caps

Magnetic Caps are made from a special polished steel alloy for optimum reliability.

Cat. No.	Description	Qty
62700	20mm Aluminum Cap, Ultra Low Bleed, Beige PTFE / White Silicone	100
62800	20mm <b>Blue Magnetic</b> Aluminum / Tin Cap, Ultra Low Bleed, Beige PTFE / White Silicone	100
62900	20mm <b>Magnetic</b> Aluminum Cap, Ultra Low Bleed, Beige PTFE / White Silicone	100

For pricing, see [analytical-sales.com](http://analytical-sales.com)



62546

62546A

62346

62346A

72546


72546A

### 18mm Screw Thread Head Space Vials

Cat. No.	Description	Qty
62546	10mL Clear Glass Flat Bottom Vial, 22.5 x 46.0mm	100
62546A	10mL Amber Glass Flat Bottom Vial, 22.5 x 46.0mm	100
62546S	10mL Clear Glass Flat Bottom Vial, 22.5 x 46.0mm - Silanized	100
62546S-CASE	10mL Clear Glass Flat Bottom Vial, 22.5 x 46.0mm - Silanized	1000
62346	10mL Clear Glass Round Bottom Vial, 22.5 x 46.0mm	100
62346A	10mL Amber Glass Round Bottom Vial, 22.5 x 46.0mm	100
62346S	10mL Clear Glass Round Bottom Vial, 22.5 x 46.0mm - Silanized	100
62346S-CASE	10mL Clear Glass Round Bottom Vial, 22.5 x 46.0mm - Silanized	1000
72546	20mL Clear Glass Flat Bottom Vial, 22.5 x 75.5mm	100
72546A	20mL Amber Glass Flat Bottom Vial, 22.5 x 75.5mm	100
72546S	20mL Clear Glass Flat Bottom Vial, 22.5 x 75.5mm - Silanized	100
72546S-CASE	20mL Clear Glass Flat Bottom Vial, 22.5 x 75.5mm - Silanized	1000


### 18mm Magnetic Screw Thread Closures with Ultra Low Bleed Silicone Liners

Cat. No.	Description	Qty
62515	18mm Magnetic Screw Cap with Teflon Liner, 0.010"	100
62518	18mm Magnetic Screw Cap with Blue Teflon / White Silicone, 0.080"	100
62517	18mm Magnetic Screw Cap with Blue Teflon / White Silicone, 0.080"	1000



### Custom Head Space Vial Holders

Cat. No.	Description	Qty
6233	Head Space Vial Holders for 18mm & 20mm Vials, 3 x 3	1
6244	Head Space Vial Holders for 18mm & 20mm Vials, 4 x 4	1
6255	Head Space Vial Holders for 18mm & 20mm Vials, 5 x 5	1
6266	Head Space Vial Holders for 18mm & 20mm Vials, 6 x 6	1



### Crimpers & Decappers

Cat. No.	Description	Qty
90200	Vial Crimper for 20mm OD Vials	1
90100	Decapping Pliers for All Vials	1

For pricing, see [analytical-sales.com](http://analytical-sales.com)