

**Economical!**

## PEEK™ Ultra High Pressure Fingertight Fittings

**MOST POPULAR**

### One-Piece Fittings




66518-3




60126X

### Two-Piece Fittings



### Three-Piece Fittings



Cat. No.	Description	Pressure Rating	Qty
66518-3	10-32 PEEK One-Piece Fitting for 1/16" OD Tubing	12,000 psi	10
60126X	6-32 PEEK One-Piece Fitting for 1/32" OD Tubing	15,000 psi	10
55195X	10-32 Short PEEK Nut, with Ferrule for 1/16" OD Tubing	9,000 psi	10
55196X	10-32 Long PEEK Nut, with Ferrule for 1/16" OD Tubing	9,000 psi	10
55192X	PEEK Ferrule for 1/16" OD Tubing	9,000 psi	10
60110X	10-32 Nut for 1/16" OD Tubing	16,500 psi	10
60100X	PEEK Ferrule for 1/16" OD Tubing w/SS Ring	16,500 psi	10
60132X	PEEK Ferrule for 1/32" OD Tubing w/SS Ring	16,500 psi	10

For pricing, see [analytical-sales.com](http://analytical-sales.com)