

**Finger Tight  
Stainless Steel Fittings  
on 1/16" OD Stainless Steel Tubing**

**Rated to 27,000 psi**



Graphite Ferrule

**Silver Fitting:** finger tighten to seat ferrule into port.

**Purple Locking Cap:** finger tighten to "lock in" tubing. Prevents tubing "creep" and creates a true ZDV connection.

**New! Stopper O-ring**

Prevents backward travel of Locking Cap when connecting or disconnecting fitting, position can be adjusted along tubing

## Ultra High Pressure Finger Tight UHPLC Assemblies

- Zero Dead Volume Connection
- Pressure Tested to 27,000 psi
- Two-part Finger Tightening System "Locks In" Ferrule and Tubing



### Replacement Graphite Ferrules

Cat. No.	Dimensions	Qty
68566-1	One-Piece Ferrule System	5



### Single UHPLC Assembly Stainless Steel Tubing with One Fitting

Cat. No.	Dimensions
61510	10cm x 0.005" ID
61710	10cm x 0.007" ID
61610	10cm x 0.010" ID
61210	10cm x 0.020" ID
61715	15cm x 0.007" ID
61520	20cm x 0.005" ID
61720	20cm x 0.007" ID
61620	20cm x 0.010" ID
61220	20cm x 0.020" ID
61525	25cm x 0.005" ID
61530	30cm x 0.005" ID
61730	30cm x 0.007" ID
61630	30cm x 0.010" ID
61230	30cm x 0.020" ID
61540	40cm x 0.005" ID
61740	40cm x 0.007" ID
61640	40cm x 0.010" ID
61240	40cm x 0.020" ID
61550	50cm x 0.005" ID
61750	50cm x 0.007" ID
61650	50cm x 0.010" ID
61250	50cm x 0.020" ID
61560	60cm x 0.005" ID
61750	70cm x 0.007" ID



### Double UHPLC Assembly Stainless Steel Tubing with Two Fittings

Cat. No.	Dimensions
62510	10cm x 0.005" ID
62710	10cm x 0.007" ID
62610	10cm x 0.010" ID
62210	10cm x 0.020" ID
62715	15cm x 0.007" ID
62520	20cm x 0.005" ID
62720	20cm x 0.007" ID
62620	20cm x 0.010" ID
62220	20cm x 0.020" ID
62530	30cm x 0.005" ID
62730	30cm x 0.007" ID
62630	30cm x 0.010" ID
62230	30cm x 0.020" ID
62540	40cm x 0.005" ID
62740	40cm x 0.007" ID
62640	40cm x 0.010" ID
62240	40cm x 0.020" ID
62550	50cm x 0.005" ID
62750	50cm x 0.007" ID
62650	50cm x 0.010" ID
62250	50cm x 0.020" ID
62560	60cm x 0.005" ID
62760	60cm x 0.007" ID
62501	100cm x 0.005" ID

**Longer lengths are available - call for pricing**

## FlexChrom<sup>®</sup>

Flexible Stainless Steel Tubing

### Assemblies

UHPLC, Pre-Finished, One-Piece Flexible Stainless Steel Tubing with Finger Tight Stainless Steel Fittings

Rated to 27,000 psi



Two-part Finger Tight Stainless Steel Fitting - "locks in" ferrule and tubing

**New! Stopper O-ring**  
Prevents backward travel of Locking Cap when connecting or disconnecting fitting, position can be adjusted along tubing

• 1/16" O.D.

• One-piece Stainless Steel – no brazed-on sleeves

• True ZDV

• Greater flexibility

• Pressure tested to 27,000 psi

• Thin, Durable Tubing Will Not Crack or Snap

One-Piece Stainless Steel

▲ 1/16" O.D.

Flexible 1/16" Stainless Steel Tubing - bends easily to fit into tight or hard to reach spaces

#### Replacement Graphite Ferrules

Cat. No.	Dimensions	Qty
68566-1	One-Piece Ferrule System	5



### 1/16" Flexible SS Tubing with SS Fingertight Fittings

#### Stainless Steel Tubing with No End Fittings

Cat. No.	Dimensions
990011-Flex	10cm x 0.005"ID
930011-Flex	10cm x 0.007"ID
970011-Flex	10cm x 0.010"ID
990015-Flex	20cm x 0.005"ID
930015-Flex	20cm x 0.007"ID
970015-Flex	20cm x 0.010"ID
990028-Flex	30cm x 0.005"ID
930028-Flex	30cm x 0.007"ID
970028-Flex	30cm x 0.010"ID



#### Single UHPLC Assembly

Cat. No.	Dimensions
<b>990141-Flex</b>	<b>10cm x 0.004"ID</b>
990111-Flex	10cm x 0.005"ID
930111-Flex	10cm x 0.007"ID
970111-Flex	10cm x 0.010"ID
<b>990142-Flex</b>	<b>20cm x 0.004"ID</b>
990115-Flex	20cm x 0.005"ID
930115-Flex	20cm x 0.007"ID
970115-Flex	20cm x 0.010"ID
<b>990143-Flex</b>	<b>30cm x 0.004"ID</b>
990128-Flex	30cm x 0.005"ID
930128-Flex	30cm x 0.007"ID
970128-Flex	30cm x 0.010"ID



#### Double UHPLC Assembly

Cat. No.	Dimensions
<b>990241-Flex</b>	<b>10cm x 0.004"ID</b>
990211-Flex	10cm x 0.005"ID
930211-Flex	10cm x 0.007"ID
970211-Flex	10cm x 0.010"ID
<b>990242-Flex</b>	<b>20cm x 0.004"ID</b>
990215-Flex	20cm x 0.005"ID
930215-Flex	20cm x 0.007"ID
970215-Flex	20cm x 0.010"ID
<b>990243-Flex</b>	<b>30cm x 0.004"ID</b>
990228-Flex	30cm x 0.005"ID
930228-Flex	30cm x 0.007"ID
970228-Flex	30cm x 0.010"ID

For pricing, see [analytical-sales.com](http://analytical-sales.com)