


## BLACK TrueTaper™

*For Light Sensitive Compounds*



*TrueTaper™ 96-Well Collection Plates were designed for chromatographers faced with limited sample quantity. TrueTaper™ incorporates a smooth, tapered well shape which allows the entire sample to gravitate to the bottom of the well.*

- Square wells with unique **TAPERED** round bottoms for maximum sample recovery
- 100uL tapered bottom volume
- 2.0mL total volume capacity

Cat. No.	Description	Qty
968820BL	TrueTaper 96-Well Collection Plate, Black	25

- Also available with **Siliguard** non-binding coating -  
*For pricing, see [analytical-sales.com](http://analytical-sales.com)*