

Canary-Safe™
Solvent Safety products

GL38 Mobile Phase Bottle Safety Caps



- Screws directly onto **GL38 threaded amber solvent bottles** (*NOTE: Most 4L solvent bottles are GL38 threads*)
- **No thread adapter needed!**
- Comes with **Air Check Valve and Fittings**

GL38 Mobile Phase Bottle Safety Caps

Cat. No.	Description
CCGL3801	Canary-Safe Cap, GL38, 1 Standard Tubing Port
CCGL3802	Canary-Safe Cap, GL38, 2 Standard Tubing Ports
CCGL3803	Canary-Safe Cap, GL38, 3 Standard Tubing Ports



Replacement Filter Check Valve

- Fits Canary-Safe™ mobile phase safety caps
- Replace every six months for optimum protection.

Cat. No.	Description
CX0001	Replacement Filter Check Valve for Mobile Phase Safety Caps