



Fingertight UPLC Fitting

Double Ferrule Stainless Steel Fitting

- **Rated to 13,500psi**
- **Unique Double Ferrule, Fingertight, ZDV**
- **The Grooved Stainless Steel Collet Reinforces the Polymer Ferrule**

| Cat. No. | Description | Qty |
|----------|--|-----|
| SDY-201 | SS Nut, includes Polymer Ferrule and SS Collet | 2 |
| SDY-102 | Replacement Polymer Ferrules | 5 |
| SDY-103 | Replacement Stainless Steel Collets | 2 |

For pricing, see analytical-sales.com