

Echelon C18: Sub 2 μ m & 3 μ m Columns



Echelon C18 provides:

- Proprietary highest quality C18 & packing method
- Better retention and resolution of compounds more polar than traditional C18 can manage
- **Extreme Durability:**
 - pH 2-11
 - Handles 100% aqueous
- **Proven Reliability:**
 - Adopted by high-throughput screening groups for years

Featured in 2013 ASMS poster, visit:

<http://www.analytical-sales.com/Echelon-Columns.html>

All this now in a UHPLC Column!

<2 μ m Echelon Columns

C18 100Å 1.7 μ m

Cat. No.	Description
ADV8106	Echelon C18, 100Å, 1.7 μ m, 20mm x 2.1mm
ADV8107	Echelon C18, 100Å, 1.7 μ m, 30mm x 2.1mm
ADV8108	Echelon C18, 100Å, 1.7 μ m, 50mm x 2.1mm
ADV8117	Echelon C18, 100Å, 1.7 μ m, 100mm x 2.1mm

3 μ m Echelon Columns

C18 100Å 3 μ m

ADV8302	Echelon C18, 100Å, 3 μ m, 20mm x 2.1mm
ADV8303	Echelon C18, 100Å, 3 μ m, 30mm x 2.1mm
ADV8304	Echelon C18, 100Å, 3 μ m, 50mm x 2.1mm

For pricing, see [analytical-sales.com](http://www.analytical-sales.com)