




• Includes fittings for 1/8" tubing



Multi-Inlet Waste Manifold


Cat. No.	Description
MW0802	Canary-Safe Multi-Inlet Waste Manifold with 8 Standard Tubing Ports for 1/16" or 1/8" OD Tubing, 2 Ports for Barbed Adapter, and 1 Large Outlet. Includes 1/8" Fittings and Barbed Adapters

Other configurations available, see analytical-sales.com




Keep your solvent waste lines neat and organized with our new Multi-Inlet Waste System Manifold!

- Made of durable chemical-resistant polypropylene
- Fits bench up to 3" thick
- Includes fittings and barbed adapters
- Other configurations available - call for details



For 1/16" or 2.5mm tubing, use these fittings (sold separately)

Cat. No.	Description	Qty
59249	PEEK™ One Piece Fitting (1/4-28) for 1/16" OD Tubing	Ea
59305X	Green Delrin Flangeless Male Nut, use with Ferrule (59353X)	10
59353X	ETFE Ferrule for 2.5mm OD Tubing	10






Plugs

Cat. No.	Description	Qty
CF5007	Plugs For 1/8" or 1/16" Ports, White PFA	10
CF5009	Plugs for 6.4mm ID Barbed Connector Port, PP	5

For pricing, see analytical-sales.com

The Canary-Safe™ Multi-inlet Waste Manifold remediates these waste solvent situations:



Exhaust Filter Adapters for JustRite® Waste Containers

Cat. No.	Description	Qty
CJX01	Exhaust Filter Fitting/ <i>Long</i> Adapter for SS Port on JustRite Waste Container	Each
CJX02	Exhaust Filter Fitting/ <i>Short</i> Adapter for Gray PEEK Port on JustRite Waste Container	Each



Canary-Safe™ Vapor Exhaust Filters for Safety Waste Caps

Cat. No.	Size	Description	Replacement	Qty
CX3000	S	Exhaust Filter, for containers < 5L	Every 3 Months	Ea
CX6000	M	Exhaust Filter, for containers < 20L	Every 6 Months	Ea
CX12000	L	Exhaust Filter, for containers > 20L	Every 12 Months	Ea
CX6010	M	Exhaust Filter with H ₂ SO ₄ (sulfuric acid) for Volatile Bases	Every 6 Months	Ea
CX6020	M	Exhaust Filter with KOH (potassium hydroxide) for Volatile Acids	Every 6 Months	Ea



Timer/Change Indicator

1. Push button to activate
2. Replace when service life expires

Exhaust Filters with Timer/Change Indicator

Canary-Safe exhaust filters contain an activated carbon, multi-component granulate material which absorbs 99% of all vapor and particulate materials evaporating from your solvents. Filters come with a **Timer/Change Indicator** at no additional cost.



Avoid Messy Overflow!
Check out our
Sentinel Fluid Level Alarm System