

Hamilton Tip Compatibility with [Select-A-Vial](#)

	<i>50μl Tip</i>	<i>300μl Tip</i>	<i>1000μl Tip</i>
<i>350μL Select-A-Vial</i>	Reaches Bottom	Reaches Bottom	Tip Cannot Reach Bottom
<i>500μL Select-A-Vial</i>	Reaches Bottom	Reaches Bottom	Reaches Bottom
<i>1mL Select-A-Vial</i>	Reaches Bottom	Reaches Bottom	Reaches Bottom