




NOTES:

- IDENTICAL TO 24625 REACTOR BLOCK ASSEMBLY BUT WITH 18mm WELL SPACING
- SCREW HOLE PATTERN ALSO CHANGED RESPECTIVE TO THE NEW TIGHTER SPACING

REVISION HISTORY			
REV	DESCRIPTION	DATE	APPROVED
1A	INITIAL RELEASE	3/18/2021	ENGINEERING
1B	ADDED DIMENSIONS	8/14/2024	ENGINEERING



DRAWN Matt Dabrowski email engineering@analytical-sales.com	3/18/2021	ANALYTICAL SALES & SERVICES, INC © 2024	
BREAK ALL EDGES		TITLE 24626_4ml_open_bot_18mm_rev1b_PUBLIC	
Tolerance (Unless Noted): .XXX = ± .005 .XX = ± .015 Angular ± 2°		SIZE B	REV 1B
		SCALE 1 : 1	SHEET 1 OF 1