


Component Schematic

96-well Microplate (PP) Deep Well (Square) 1ml, 2ml & 350uL



| | | | |
|----|---------------|------|--------------|
| L | 127.60 ± 0.10 | H350 | 15.30 ± 0.10 |
| W | 85.50 ± 0.10 | H1 | 24.50 ± 0.10 |
| H | 44.40 ± 0.05 | H2 | 39.90 ± 0.05 |
| Hs | 4.00 ± 0.10 | P | 9.00 ± 0.05 |
| X1 | 14.00 ± 0.05 | D | 8.25 ± 0.05 |
| Xn | 113.50 ± 0.05 | C | 12 |
| Y1 | 11.00 ± 0.05 | R | 8 |
| Yn | 74.50 ± 0.05 | | |

Volume 350µl = 750µl, 1ml = 1,400µl, 2ml = 2,500µl
 Well form V Bottom